

Light Duty Cautery Sets-ENT



Features

- Modern 230 volt mains transformer designed for optimum energy efficiency
- Compact and lightweight
- Completely portable with soft carry case & shoulder strap - unit can also be used whilst in the case
- Unique slot on front of unit safely stores handle and burner when not in use (for use with Pen Handle JA123 only)
- Easy to use with illuminated ON/OFF button
- Extensive range of reusable and disposable burners available in 3 lengths
- Disposable handle covers available which reduce the risk of infection.

JB104 Cautery Handle Cover



JC620

- LD Cautery Set consisting of JA620 Transformer, JA123 LD Cautery Handle and 5 Burners JA321D, JA321G, JA321I, JA321O and JA321R.
- This cautery set is designed for use in confined spaces. The cautery burners are suitable for use in the wet tissue areas of the nose and mouth.



CAUTERY