

Light Duty Cautery Set - Ophthalmic



Features

- Modern 230 volt mains transformer designed for optimum energy efficiency
- Compact and lightweight
- Completely portable with soft carry case & shoulder strap - unit can also be used whilst in the case
- Easy to use with illuminated ON/OFF button
- Range of reusable and disposable burners
- Disposable handle covers available which reduce the risk of infection

JC623

Ophthalmic Cautery Set

230 volt LD Ophthalmic Cautery set consisting of Transformer, JA114 Eye Handle and Eye Burners JA211A, JA211B, JA211F, JA211S and JA211CP



Cautery