

Light Duty Cautery Sets - Dermatology



Features

- For dermatology related minor operations
- Modern 230 volt mains transformer designed for optimum energy efficiency
- Compact and lightweight
- Completely portable with soft carry case & shoulder strap - unit can also be used whilst in the case
- Unique slot on front of unit safely stores handle and burner when not in use (for use with Pen Handle JA123 only)
- Easy to use with illuminated ON/OFF button
- Complete with 5 reusable cautery burners: coagulation ball (JA321B), ring cutter (JA321C), blade cutter (JA321F), standard point (JA321S) and cold point (JA321CP)
- 2 other handle types available: Pistol Grip (JA143) and Mark Hovell (JA101)
- Extensive range of reusable and disposable burners available in 3 lengths
- Disposable handle covers available which reduce the risk of infection.

JC620	LD Cautery Set consisting of JA620 Transformer, JA123 LD Cautery Handle and 5 Burners JA321B, JA321C, JA321F, JA321S and JA321CP.
JC620-1	LD Cautery Set consisting of JA620 Transformer, JA123 LD Cautery Handle & 20 S/Use burners
JC621	LD Cautery Set consisting of JA620 Transformer, JA123 LD Cautery Handle – No Burners
JC622	LD Cautery Set consisting of JA620 Transformer, JA101 Mark Hovell Handle, JB123 LD Cautery Cord with shrouded pins, JA321B, JA321C, JA321F, JA321S and JA321CP.
JC623	LD Ophthalmic Cautery Set consisting of JA620 Transformer, JA114 Eye Handle and Eye Burners JA211A, JA211B, JA211F, JA211S and JA211CP
JC624	LD Cautery Set consisting of JA620 Transformer JA143 LD Pistol Grip Handle and 5 Burners JA321B, JA321C, JA321F, JA321S and JA321CP.
JC625	LD Cautery Set consisting of Transformer, JA143 Pistol Grip Handle, no burners
JC626	LD Cautery Set consisting of Transformer, JA143 LD Pistol Grip Handle and 5 LD Burners of your choice



Cautery