

Nasopharyngoscopes

SW3012



Flexible Naso-Pharyngoscope with conventional image bundle, 7000 Pixel, 2-X Tip Deflection, One Hand operation, Designed for gas sterilisation.

Working Length:	300 mm
Working \varnothing :	3.2 mm
Tip	Deflection: Up 125°, Down 125°
Tip Section:	20mm
Viewing Direction:	0° Wide Angle

SW3011



Flexible Naso-Pharyngoscope with conventional image bundle, 12000 Pixel, 2-X Tip Deflection, One Hand operation, Designed for gas sterilisation.

Working Length:	300 mm
Working \varnothing :	3.8 mm
Tip	Deflection: Up 160°, Down 160°
Tip Section:	25mm
Viewing Direction:	0° Wide Angle