

Trans-rectal Guided Injection Needle

Order Code: SERA2104

1.2mm Ø x 250mm long (min 10)

Biopsies are commonly performed from ultrasound guidance. During the procedure, a special biopsy needle is inserted into the prostate gland through the wall of the rectum to remove small samples of tissue for pathologic analysis. This method is known as trans-rectal ultrasound (TRUS) guided biopsy.

RB Medical can supply a disposable injection needle for use in TRUS guided procedures to diagnose prostate cancer using a Trans-rectal scanner.



Urology