

# LLETZ Specula

- Designed for 'LLETZ' Procedures (Large Loop Excision of the Transformation Zone)
- Teflon coated to ensure highest levels of protection for patients and physicians
- Used with electrosurgical units
- Provision to connect with a smoke evacuator unit

## 4-Way Lateral Expanding LLETZ

Combines the characteristics of a standard speculum and a vaginal wall retractor



DA211	4-way Cusco speculum, Teflon coated	Large
DA212	4-way Cusco speculum, Teflon coated	Medium

## Standard LLETZ



DA130-107	Graves speculum, smoke flume, facing up	Extra large
DA133-107	Pederson speculum, 25mm opening x 170mm blade left side open	Extra large
DA133-108	Pederson speculum, 25mm x 140mm blade, left side open	Extra large

Cont....

# LLETZ Specula

## LLETZ Vu-More



DA132-104	Graves Vu-More, 40mm opening x 150mm, open tube	Extra large
DA132- 107	Graves Vu-More, 40mm opening x 180mm, open tube	Extra large
DA132-1071	Graves Vu-More, 40mm opening x 180mm, closed tube	Extra large

We can provide other types of LLETZ instruments, please discuss with your local Area Development manager.