

Operating and Diagnostic Hysteroscopy

Solutions from RB Medical

Sheaths

- Sheaths to fit most makes of hysteroscope
- Reusable and single use versions
- Reliable irrigation
- Single and continuous flow
- Working channel permits use of 5fr instruments and bipolar electrodes
- Sheath enables complete rotation of scope for convenient viewing



RA229C – for use with 2mm Hysteroscope

- Twin Channel
- Ratchet Handle on Grasper
- Non-Ratchet on Scissor
- Extended working length



Hysteroscopes

- Where diameter is critical for safe inspection of the uterine cavity, we offer
- 4 sizes of hysteroscope:
 - 2 mm – 0 degree
 - 2.9 mm – 0, 12 and 30 degree
 - 4.0 mm – 0, 12, 30 and 70 degree
 - 5.0 mm – 0, 30 degree
- Small diameter still offers excellent image quality
- Comprehensive range of sheaths
- Compatible with most cameras and light sources
- Fully autoclavable with outstanding record of maintained image quality
- Protection tubes & sterilising baskets/cages available

Instruments

Full range of biopsy cups, grasping forceps and scissors available



Gynaecology