

Spackmann Cannula

SE5054 Spackmann Cannula 30cm Tube Length (Min 20)

100% infection control by using single use

Same shape and size as reusable instrument

Large diameter stem for improved handling and ease of use

Theatre quality stainless metal instrument

- Fixed anatomical cervical cone provides leak free seal
- Atraumatic tip provides reduced uterine perforation risk
- Even distribution of injected dyes and fluids from 13 holes on tip



Adjustable clamp with tenaculum attachment

Handle has tube and Luer taper fittings for injection of dye into the uterine cavity