

NOW ONLY AVAILABLE FROM RBMEDICAL

Order Code YA501

Special Offer

Buy 2 boxes at a time £100 (per Box 25) and FREE delivery + add Loops & Balls, or any other instruments to your order and get them delivered FOC too!

A unique, sterile, disposable curette for sampling the uterine endometrium

- **Only one EndoSampler™ required per patient** – vacuum can be recreated whilst leaving curette in place
- **Easy transfer of sample** -backflow is prevented by lock spring mechanism on syringe
- **Prevents contamination during withdrawal** - sample is collected and secured prior to curette removal
- **Better sampling** – flexible tube is conveniently bent bringing it immediately into contact with the wall of the uterus



Graduated markings on cannula can be easily read



Attaching syringe to the base of curette

Semi-rigid tapered 3mm cannula enables easier introduction and minimises the possibility of perforation

Facilitates both **endometrial scraping and aspiration** for comprehensive specimen removal at the same time.