

MedGyn Pipette IV (Four Holes) Endometrial Sampler YA505

MedGyn Pipette IV is used for obtaining endometrial tissue in the outpatient setting. It is designed to improve ease of tissue collection and to reduce the physical discomfort to the patient.

- Polypropylene sheath 3.1 mm external diameter
- 2.6mm internal diameter
- 23.5cm long
- 4 lateral orifices - 2.5mm in diameter
- 11 guide marks



Features

- Sterile, single use, latex free
- Semi-rigid, tapered 3.1mm cannula enables easier insertion and minimises the possibility of perforation
- Four holes in distal end for improved sample collection of the uterine cavity
- Easier dispensing of histology specimen – plunger slides freely inside cannula and is easy to push even when fully extended
- Flexible cannula offers better sampling from all areas of the endometrial cavity including difficult to reach cornual areas
- Can be manipulated to retain a fixed curved shape
- Stop on proximal end of cannula prevents plunger overrunning and jamming at distal end
- Ergonomically designed thumb plate