

LED Light Source



GA472

- Portable and compact
- 50,000 hours life LED
- Cold, white, high intensity light
- Universal fibre optic connection
- Electronic dimmer

Technical specifications

LED	80w
Light colour temperature	6000°k
Power consumption	110 watts
Class	1
Type	BF
LED life	50,000 hours approximately
Over temperature protection	Automatic LED system power-off
Dimmer	Electronic
Fibre optic connector	Universal, compatible with major brands
Dimension	202 x 290 x 95mm (W x D x H)
Weight	2.5kg
Power Supply	100-240V
Frequency	50/60Hz



Light Sources