

RB Medical offer endoscopes with fittings compatible with all major manufacturers locking collars. Our rigid endoscopes have the highest quality optical features and are built with the following technical benefits:

- Distal and proximal sapphire windows to reduce washer corrosion
- Water absorbing filters to reduce hot / cold condensation
- Rod lenses in short triplets to reduce fragile nature of scope
- Thermal compensating springs to reduce autoclave stress
- Temperature sensors to check for over heating
- Light stem red heat filter and light amplifier
- Stainless steel outer case, laser welded to give uniform expansion



**SW3002** 0 degree Neuroscope 2.7mm $\varnothing$  , 213 mm Length

**SW3003** 30 degree Neuroscope 2.7mm  $\varnothing$  x, 213 mm Length

**SW2126** 0 Degree Neuroscope 2.7mm  $\varnothing$  x187mm Length

**SW2130** 0 Degree Neuroscope 2.7mm  $\varnothing$  x 187mm Length