


Lumbar Spine procedures

McCulloch Retractor Set

Reduce Cross-Infection



McCulloch Retractor Set

SE559135 McCulloch Retractor Set

- McCulloch retractor set for use in lumbar spine procedures with disposable body and accessories
- Blades slide smoothly onto the frame for easy set up and break down
- A variety of blades offered including toothed, hook and muscle in a variety of sizes/widths

Available as a complete set, or, select a set of blades of your choice from the following accessories:



Neurosurgery

Lumbar Spine procedures

Muscle Blades

SE559135-50N Narrow muscle blade 50mm



SE559135-50N-S Narrow muscle blade set (x2) 50mm



SE559135-50W Wide muscle blade 50mm



SE559135-50W-S Wide muscle blade set (x2) 50mm



SE559135-60N Narrow muscle blade 60mm



SE559135-60N-S Narrow muscle blade set (x2) 60mm



Lumbar Spine procedures

SE559135-60W Wide muscle blade 60mm 	SE559135-60W-S Wide muscle blade set (x2) 6.0mm 
SE559135-70N Narrow muscle blade 70mm 	SE559135-70N-S Narrow muscle blade set (x 2) 70mm 
SE559135-70W Wide muscle blade 70mm 	SE559135-70W-S Wide muscle blade set (x 2) 70mm 
SE559135-80N Narrow muscle blade 80mm 	SE559135-80N-S Narrow muscle blade set (x 2) 80mm 

Lumbar Spine procedures

SE559135-80W Wide muscle blade 80mm



SE559135-80W-S Wide muscle blade set (x 2) 80mm



Hook Blades

SE559135-HE Single hook blade Size 3



SE559135-H4 Single hook blade Size 4



SE559135-H5 Single hook blade Size 5



SE559135-H6 Single hook blade Size 6



Toothed Blades

SE559135-H-S Single hook blade Set



SE559135-M3 Multiple toothed blade Size 3



SE559135-M4 Multiple toothed blade Size 4



SE559135-M5 Multiple toothed blade Size 5



Lumbar Spine procedures

SE559135-M6 Multiple toothed blade Size 6



SE559135-M-S Multiple toothed blade Set



Available as a complete set, or select a set of your choice from the accessories above: