


# Biopsy & Diagnostic



## Operating Endoscopes & Sheaths

| Endoscope  | Sheath                           |   | W/L              | Size         | Inst. Size     |
|--|----------------------------------|---|------------------|--------------|----------------|
| SW4004 0° Scope Dia 4mm Length 302mm – Wide Angle  | RC235-0<br>RA257-0<br>RA258-0    | Double Channel - Laser<br>Inner Single Channel<br>Outer | 200mm<br>201.5mm | 26Fr<br>17Fr | 2 x 5Fr<br>5Fr |
| SW4005 12° Scope Dia 4mm Length 302mm – Wide Angle | RC235-12<br>RA257-12<br>RA258-12 | Double Channel - Laser<br>Inner Single Channel<br>Outer | 200mm<br>201.5mm | 26Fr<br>17Fr | 2 x 5Fr<br>5Fr |
| SW4012 30° Scope Dia 4mm Length 302mm              | RC235<br>RA257<br>RA258          | Double Channel - Laser<br>Inner Single Channel<br>Outer | 200mm<br>201.5mm | 26Fr<br>17Fr | 2 x 5Fr<br>5Fr |

## Biopsy Endoscopes & Sheaths

| Endoscope  | Sheath   | W/L   | Size |
|--|----------|-------|------|
| SW4004 0° Scope Dia 4mm Length 302mm – Wide Angle  | RA234-0  | 261mm | 5mm  |
| SW4005 12° Scope Dia 4mm Length 302mm – Wide Angle | RA234-12 | 261mm | 5mm  |
| SW4012 30° Scope Dia 4mm Length 302mm              | RA234    | 261mm | 5mm  |
| SW4016 70° Scope Dia 4mm Length 302mm              | RA234-70 | 261mm | 5mm  |

### 5Fr Instruments:

|         |  |         |  |
|---------|--|---------|--|
| SW5005  | 5Fr Flexible Biopsy Forceps Double Action Spoon Jaws | SW5023  | 5Fr Flexible Grasping Forceps – Un-insulated |
| SW5006  | 5 Fr Biopsy Forceps length 600mm                     | SW5023R | 5Fr Flexible Grasping Forceps – Semi Rigid   |
| SW5019  | 5Fr Biopsy Scissors                                  | SW5024  | 5Fr Grasping Forceps                         |
| SW5019R | 5Fr Flexible Scissors-semi rigid                     | SW51005 | 5Fr Coagulation Button Electrode             |
| SW5020  | 5Fr Flexible Scissors – Semi Rigid                   | SW51015 | 5Fr Needle Electrode                         |

### 3Fr Instruments:

|        |                                 |         |                               |
|--------|---------------------------------|---------|-------------------------------|
| SW4074 | 3Fr Flexible Scissors           | SW4082  | 3Fr Flexible Grasping Forceps |
| SW4076 | 3Fr Flexible Grasping Forceps   | SW4048  | 3Fr Double Biopsy Forceps     |
| SW4078 | 3Fr Double Spoon Biopsy Forceps | SW51003 | 3Fr HF Button Electrode       |



## Urology