


Cystoscopy/Urethroscopy



Endoscopes		Sheaths
SW4004 0° Scope Dia 4mm Length 302mm – Wide Angle	SX4040	16Fr Sheath including obturator
SW4005 12° Scope Dia 4mm Length 302mm – Wide Angle	SX4046	19Fr Sheath including obturator
SW4012 30° Scope Dia 4mm Length 302mm	SX4020	22Fr Sheath including obturator
SW4016 70° Scope Dia 4mm Length 302mm	SX4052	24Fr Sheath including obturator

Catheter Deflecting Mechanisms

- SW4026 Catheter deflecting mechanism including Albarran lever, 1 working channel to fit sheaths, 19,22 and 24Fr
 SW4028 Catheter deflecting mechanism including Albarran lever, 2 working channels to fit sheaths, 19,22 and 24Fr

Catheter Capacities

Sheaths with Obturators	Deflecting SW4026	Mechanism SW4028	Telescope Bridge
SX4040 16Fr	-	-	1 x 5Fr
SX4046 19Fr	1 x 7Fr	2 x 5Fr	1 x 7Fr
SX4020 22Fr	1 x 7Fr	2 x 5Fr	1 x 7Fr
SX4052 24Fr	1 x 7Fr	2 x 5Fr	1 x 7Fr

Catheter Introducers

- SW4094 Guyon catheter Introducer
 SW4094L Guyon catheter introducer (long)
 SW4094L-1 Guyon catheter introducer (long) Macclesfield Design
 SW4094L-2 Guyon catheter introducer (long) Frimley Design
 SW409L-4 Guyon catheter introducer (extra - long) Frimley Design
 SW4094-1 Guyon catheter introducer

Cystoscopy/Urethroscopy

Cystoscope Needles

SERA2101 Disposable 19g needle

FLEXIBLE INSTRUMENTS

5Fr Instruments:

SW5005	5Fr Flexible Biopsy Forceps Double Action Spoon Jaws	SW5023	5Fr Flexible Grasping Forceps – Un-insulated
SW5006	5Fr Biopsy Forceps Length 600mm	SW5023R	5Fr Flexible Grasping Forceps – Semi Rigid
SW5019	5Fr Flexible Scissors	SW5023L	5Fr flexible grasping forceps – semi rigid 600mm length
SW5019R	5Fr Flexible Scissors – Semi Rigid	SW5023RL	5Fr flexible grasping forceps – semi rigid 600mm length
SW5020	5Fr Scissors Length 450mm	SW5024	5Fr Grasping Forceps Length 450mm
			5Fr Coagulation Button Electrode
		SW51005	5Fr Needle Electrode

3Fr Instruments:

SW4074	3Fr Flexible Scissors Length 510mm	SW4082	3Fr Flexible Grasping Forceps Length
SW4076	3Fr Flexible Grasping Forceps Length 510mm	SW4048	3Fr Double Biopsy Forceps Length
SW4078	3Fr Double Spoon Biopsy Forceps Length 510mm	SW51003	3fr HF Button Electrode



Urology