

FlexiLine scopes are low cost system built telescopes that do not compromise the quality of the product and the optical and sealing quality is equal to the RB Medical standard Endoscopes. Optical modules are pre-focused and sealed in highly monitored cleanroom conditions. The finish and sealing of all RB Medical endoscopes ensures the integrity of the product.

What you will get from a new *FlexiLine* endoscope

- Perfectly balanced lens system
- Factory focused to give best picture quality
- Robust construction
- Compatibility with many instrument systems



SW4004-1	<i>FlexiLine</i> 0° Scope Ø4mm Length 302mm - Wide Angle
SW4005-1	<i>FlexiLine</i> 12° Scope Ø4mm Length 302mm – Wide Angle
SW4060-1	<i>FlexiLine</i> 25° Scope Ø4mm Length 298mm – Wide Angle – Wolf fitting
SW4012-1	<i>FlexiLine</i> 30° Scope Ø4mm Length 302mm
SW4016-1	<i>FlexiLine</i> 70° Scope Ø4mm Length 302mm
SW4062-1	<i>FlexiLine</i> 70° Scope Ø4.0mm x 303mm Wolf fitting