


Long resectoscope set



Meet the needs of larger patients with the RB Medical long resectoscope set

Available for purchase or hire

Benefits:

- Tried and tested equipment
- Improved management of case lists and patient care
- No need for capital outlay when budgets are tight

Set contains:

Product Code	Description
SW4012L	30 degree telescope 4mm x 367mm long, wide angle as standard
SX5003-1L	26 Fr continuous flow sheath: inner and outer sheaths, taps with luer fittings, ceramic tip to inner sheath – rotating, standard solid obturator
SW4034L	Active working element with 2mm cable connection, standard cable to 8mm valley lab big pin
RA250L	Visual obturator with 5 Fr instrument channel for SW4012L (Schmidt type)
375-045	Monopolar cable for resectoscopes (Valleylab/Conmed) 8mm pin
GA134	Storz to Storz light cable
SX11649-63	Basket to hold instruments for washing and sterilizing

Reusable	Sterile single use	Description
X-RB1GL	RB1GL	24 Fr cutting loop electrode - long
X-RB1NL	RB1NL	24 Fr ball electrodes - long
X-RB1LL	RB1LL	24 Fr Collins Knife electrodes - long



Urology