


Optical Biopsy Forceps

- Ideal for foreign body removal and biopsy procedures
- A telescope joined with forceps provides a close-up view of the object being grasped
- The optical forceps can be used in many types of airway management procedures where accurate depth perception is imperative for visualisation



SW40401	Biopsy Forceps, double action spoon shaped	0°, 12° (30°)
SW40406	Biopsy Forceps, through cutting	0°, 12° (30°)
SW40403	Grasping Forceps, double action	0°, 12° (30°)
SW40405	Optical Scissors, double action	0°, 12° (30°)
SW40402	Biopsy Forceps, double action spoon shaped – short jaw	12°, 30°
SW404021	Biopsy Forceps, double action spoon shaped – short jaw	70°
SW40404	Optical Stone Crushing Forceps	0°, 12° (30°)