

Camera Stack System

Suitable for Theatre & Clinics

Order code SW0207-2 (with LED L/S)



Trolley with 1 drawer and 1 shelf
With transformer

Endodigi

The Endodigi is the complete digital documentation system for videos and images.

Simple and secure web sharing option. The data can also be shared through the virtual network, and accessed from any remote location.

Endodigi Camera System

Can be used with all scopes, rigid & flexible.

Light Source options

State of the art high performance LED light source. Superior Light output, Safety temperature automatic control technology. (SW0207-2)

Colour Video Printer

The UP-DR80MD is a full page medical grade digital printer designed specifically for medical colour printing applications where outstanding image quality and accurate colour reproduction are essential. Attaching directly to endoscopic or ultrasound systems via USB 2.0, this fast and reliable dye sublimation printer is compact enough for cart integration.

With the optional network adapter kit, it is robust enough to serve as a network printer for multiple rooms. The combination of the 301dpi thermal head and the premium media enable the UP-DR80MD to print images in full HD resolution with outstanding photo-realistic image quality. Prints from the UP-DR80MD are also highly durable, as the media includes built-in protective lamination.



Gynaecology