


# Diathermy Colposcopy Loops, Needle & Ball Electrodes – Sterile, Single Use

## Conisation Electrodes



Description	Width (mm)	Depth (mm)	Length (mm)	Order Code	Quantity
Diathermy Conisation Electrode Small	8	15	130	SE150-195-08	Box 20
Diathermy Conisation Electrode Medium	11	15	130	SE150-195-11	Box 20
Diathermy Conisation Electrode Medium Extended	11	21	130	SE150-198-11	Box 20

## Y Shaped Electrodes




Description	Width (mm)	Depth (mm)	Length (mm)	Order Code	Quantity
Diathermy Loop Ø 10x10 Y Shaped	10	10	130	SE150-175-10Y	Box 20
Diathermy Loop Ø 15x10 Y Shaped	15	10	130	SE150-175-15Y	Box 20
Diathermy Loop Ø 12x15 Y Shaped	12	15	130	SE150-177-15Y	Box 20
Diathermy Loop Ø 15x12 Y Shaped	15	12	130	SE150-180-12Y	Box 20
Diathermy Loop Ø 15x15 Y Shaped	15	15	130	SE150-180-15Y	Box 20
Diathermy Loop Ø 18x15 Y Shaped	18	15	130	SE150-185-15Y	Box 20
Diathermy Loop Ø 20x10 Y Shaped	20	10	130	SE150-190-10Y	Box 20
Diathermy Loop Ø 20x15 Y Shaped	20	15	130	SE150-190-15Y	Box 20
Diathermy Loop Ø 20x18 Y Shaped	20	18	130	SE150-190-18Y	Box 20
Diathermy Loop Ø 20x20 Y Shaped	20	20	130	SE150-190-20Y	Box 20
Diathermy Loop Ø 20x12 Y Shaped	20	12	130	SE150-190-12Y	Box 20

# Diathermy Colposcopy Loops, Needle & Ball Electrodes – Sterile, Single Use



Order Code	Description	Quantity
SE150-100	Disposable Finger Switch	Min 50
SE150-101	Disposable Finger Switch—with extended blade	Min 50
SE150-115	Diathermy Blade 3 x 17 x 130mm	Box 50
SE150-115(100)	Diathermy Blade 3 x 17 x 130mm	Box 100
SE150-120	Diathermy Needle 1 x 70mm	Box 50
SE150-125	Diathermy Needle 1 x 130mm	Box 50
SE150-120-06	Cotswold Needle 10 x 55mm, Straight	Min 20
SE150-125-01	Diathermy Needle 1 x 130mm	Box 50
SE150-125-06	Cotswold Needle 5 x 55m 45° angle	Min 20
SE150-170-05	Diathermy Loop—Square 10W x 5D x 130mm	Box 50
SE150-170-07	Diathermy Loop—Square 10W x 7D x 130mm	Box 50
SE150-170-10	Diathermy Loop-Square 10W x 10D x 130mm	Box 50
SE150-125-01(40)	Diathermy needle angled 1 x 103mm	Box 50
SE150-130	Diathermy ball electrode 2mm x 70mm	Box 50
SE150-140	Diathermy ball electrode 3mm x 60mm	Box 50
SE150-145	Diathermy ball electrode 3mm x 130mm	Box 50
SE150-150	Diathermy ball electrode 4mm x 60mm	Box 50
SE150-155	Diathermy ball electrode 4mm x 130mm	Box 50
SE150-160	Diathermy ball electrode 5mm x 60mm	Box 50
SE150-165	Diathermy ball electrode 5mm x 130mm	Box 50
SE150-175-10	Diathermy Loop 10W x 10D x 130mm L	Box 50
SE150-175-15	Diathermy Loop 10W x 15D x 130mm	Box 50
SE150-180-12	Diathermy Loop 15W x 12D x 130mmL	Box 50
SE150-180-15	Diathermy Loop 15W x 15D x 130mmL	Box 50
SE150-190-13	Diathermy Loop 20W x 13D x 130mmL	Box 20
SE150-190-15	Diathermy Loop 20W x 15D x 130mmL	Box 50
SE150-190-20	Diathermy Loop 20W x 20D x 130mmL	Box 50
SE150-200-15	Diathermy Loop 25W x 15D x 130mm	Box 20
SE150-200-20	Diathermy Loop 25W x 20D x 130mm	Box 20

# Diathermy Colposcopy Loops, Needle & Ball Electrodes – Sterile, Single Use



## Monopolar Finger Switches

Order Code	Description	Quantity
SE150-100	Single use disposable Finger Switch with blade as illustrated	Min 50
SE150-100(50)	Single Use disposable Finger Switch with blade as illustrated	Box 50
SE150-101	Single use disposable Finger Switch with extended blade	Min 50

Certified according to 93/42/EEC & compatible with generator types  
Erbe, ACC/ICC+V10, Bowa, Conmed & Valleylab