

MedGyn Silicone Pessaries

Pessaries are an effective tool for the non-surgical treatment of a variety of gynecologic problems. Commonly used in management of pelvic support defects, pessaries are also used to treat stress urinary incontinence. MedGyn offers an extensive line of quality silicone pessaries and an easy to use ring and shaatz pessary fitting set.



Ring without Support
Pessary Fitting Set



Ring with Support
Pessary Fitting Set



Shaatz Pessary
Fitting Set

The Pessary Fitting Set is an ideal diagnostic tool to help the clinician determine the proper style and size of pessary for each patient. Stored in an attractive case, the set includes six popular Ring without support, Ring with support and Shaatz pessaries and an embossed cross-reference chart to use in selecting the appropriate pessary. One set may be used for multiple patients as the latex-free rings are autoclaveable.



The Ring Pessary, available with or without a supporting membrane, is primarily used for mild first-degree uterine prolapse or a cystocele.



The Knobbed Ring pessary, available with or without a support membrane, has the same indications as the Ring pessary and may also relieve stress urinary incontinence.



The Modified Cup Pessary is primarily used for mild first-degree uterine prolapse or a cystocele.



The Gellhorn Pessary with drain feature is used to relieve the symptoms of a cystocele or a rectocele as well as a second or third-degree uterine prolapse. The pessary stem with knob is easily folded for insertion and removal.



The Gellhorn Short Stem Pessary with drain feature is used to relieve the symptoms of a cystocele or a rectocele as well as a second or third-degree uterine prolapse. The Short-Stem may provide a better fit option for patients.



The Donut Pessary may be used to relieve the stress of a second or third-degree uterine prolapse, a cystocele or rectocele. This pessary, often used by post-menopausal women, is easily compressed for insertion and removal.



The Shaatz Pessary with drain may relieve the symptoms of a cystocele as well as a second or third-degree uterine prolapse. Its convex shape ensures a snug fit while the drainage holes permit easy removal.



The Cube Pessary, available with or without drains, may be used to relieve the symptoms of a cystocele or a rectocele as well as a second or third-degree uterine prolapse.



The Cup Pessary, available with a support membrane, is used for a first or a mild second-degree uterine prolapse.



The Dish pessary, available with a supporting membrane, is used to relieve the stress of a second or third-degree uterine prolapse, a cystocele or rectocele. This pessary may also relieve stress urinary incontinence.



The Oval Pessary, available with a supporting membrane, is commonly used for a mild first-degree uterine prolapse or a cystocele. This pessary's shape inhibits its rotation in women with a history of vaginal surgery.



The Marland Pessary is available with or without support membrane, is used to relieve the symptoms of a second or third-degree uterine prolapse, a cystocele or rectocele.



The Gehrung Pessary is available with support membrane is primarily used for a significant cystocele or rectocele and may also provide relief for second or third-degree uterine prolapse.



The Hodge Pessary is available with or without support membrane, commonly used for a first or mild second-degree uterine prolapse and may also relieve a cystocele. A Hodge Pessary with support is also effective in the prevention of exercise incontinence.

Pessary	1st - 2nd Degree	2nd - 3rd Degree	Sexually Active	Vaginal Wall	Stress Urinary Incontinence	Cystocele	Rectocele	Incompetent Cervix
Cube		X	X	X	X	X	Mild X	
Cup w/ Support	X							
Modified Cup		X				X		X
Donut		X				Mild X		
Dish w/ Support	X				X	Mild X		
Gellhorn		X				Mild X		
Oval w/ Support	X					X		
Ring	X		X					
Ring w/ Support	X		X			Mild X		
Ring w/ Knobbed	X		X		X	Mild X		
Ring w/ Knobbed & Support	X		X		X	Mild X		
Shaatz w/ Drain	X					Mild X		
Marland w/ Support		X				X	X	
Marland w/o Support		X				X	X	
Hodge w/ Support	X					X		
Hodge w/o Support	X					X		

Pessary sets

<u>SKU</u>	<u>Description</u>	<u>SC NPC</u>
051000	Ring w/o Support Pessary Fitting Set of 6	
051005	Ring w/ Support Pessary Fitting Set of 6	
051006	Shaatz Pessary Fitting Set of 6	

Ring

<u>SKU</u>	<u>Description</u>	<u>SC NPC</u>
050026	Ring w/ support #1, 51 mm, 2 in	FVC1721
050027	Ring w/ support #2, 57 mm, 2.2 in	FVC1744
050028	Ring w/ support #3, 62 mm, 2.4 in	FVC1790
050029	Ring w/ support #4, 70 mm, 2.75 in	FVC1761
050030	Ring w/ support #5, 74 mm, 2.9 in	FVC1784
050031	Ring w/ support #6, 82, mm, 3.2 in	FVC1709
050032	Ring w/ support #7, 89 mm, 3.5 in	FVC1646
050033	Ring w/ support #8, 96 mm, 3.8 in	FVC1770
050034	Ring w/ support #9, 102 mm, 4 in	FVC1601
050024	Ring w/o support #1, 51mm, 2 in	FVC1694
050017	Ring w/o support #2, 57 mm, 2.2 in	FVC1627
050018	Ring w/o support #3, 62 mm, 2.4 in	FVC1780
050019	Ring w/o support #4, 70 mm, 2.75 in	FVC1725
050020	Ring w/o support #5, 74 mm, 2.9 in	FVC1661
050021	Ring w/o support #6, 82 mm, 3.2 in	FVC1811
050022	Ring w/o support #7, 89 mm, 3.5 in	FVC1613
050023	Ring w/o support #8, 96 mm, 3.8 in	FVC1600
050025	Ring w/o support #9, 102 mm, 4 in	FVC1638

Ring with knob

<u>SKU</u>	<u>Description</u>	<u>SC NPC</u>
050026K	Ring w/ knob - w/ support #1, 51mm, 2 in	FVC1732
050027K	Ring w/ knob - w/ support #2, 57mm, 2.2 in	FVC1793
050028K	Ring w/ knob - w/ support #3, 63mm, 2.5 in	FVC1738
050029K	Ring w/ knob - w/ support #4, 70mm, 2.8 in	FVC1798
050030K	Ring w/ knob - w/ support #5, 76mm, 3 in	FVC1701
050031K	Ring w/ knob - w/ support #6, 83mm , 3.3 in	FVC1711
050032K	Ring w/ knob - w/ support #7, 89mm,3.5 in	FVC1599
050017K	Ring w/ knob - w/o support #1, 51 mm, 2 in	FVC1809
050018K	Ring w/ knob - w/o support #2, 57mm, 2.2 in	FVC1624
050019K	Ring w/ knob - w/o support #3, 63mm, 2.5 in	FVC1758
050020K	Ring w/ knob - w/o support #4, 70mm, 2.8 in	FVC1608
050021K	Ring w/ knob - w/o support #5, 76mm, 3 in	FVC1772
050022K	Ring w/ knob - w/o support #6, 83mm, 3.3 in	FVC1734
050023K	Ring w/ knob - w/o support #7, 89mm,3.5 in	FVC1656

Modified cup/Cerclage

<u>SKU</u>	<u>Description</u>	<u>SC NPC</u>
050210M	Modified Cup w/ drain #1, 65mm, 2.6 in	FVC1777
050211M	Modified Cup w/ drain #2, 70mm, 2.8 in	FVC1663
050212M	Modified Cup w/ drain #3, 65mm, 2.6 in	FVC1739
050213M	Modified Cup w/ drain #4, 70mm, 2.8in	FVC1651
050214M	Modified Cup w/ drain #5, 65mm, 2.6 in	FVC1643
050215M	Modified Cup w/ drain #6, 70mm, 2.8 in	FVC1702
050191M	Modified Cup w/o drain #1, 65mm, 2.6in	FVC1757
050192M	Modified Cup w/o drain #2, 70mm, 2.8in	FVC1802
050193M	Modified Cup w/o drain #3, 65mm, 2.6 in	FVC1705
050194M	Modified Cup w/o drain #4, 70mm, 2.8 in	FVC1645
050195M	Modified Cup w/o drain #5, 65mm, 2.6 in	FVC1737
050196M	Modified Cup w/o drain #6, 70mm, 2.8 in	FVC1807

Gellhorn

<u>SKU</u>	<u>Description</u>	<u>SC NPC</u>
050155	Gellhorn w/ drain #1, 44mm, 1.7 in	FVC1728
050156	Gellhorn w/ drain #2, 51mm, 2 in	FVC1763
050157	Gellhorn w/ drain #3, 57mm, 2.2 in	FVC1733
050158	Gellhorn w/ drain #4, 63mm, 2.5 in	FVC1683
050159	Gellhorn w/ drain #5, 70mm, 2.8 in	FVC1776
050160	Gellhorn w/ drain #6, 76mm, 3 in	FVC1689
050161	Gellhorn w/ drain #7, 82mm, 3.2 in	FVC1796
050162	Gellhorn w/ drain #8, 89mm, 3.5 in	FVC1598

Gellhorn short stem

<u>SKU</u>	<u>Description</u>	<u>SC NPC</u>
050155S	Gellhorn Short Stem #1, 44mm, 1.7 in	FVC1650
050156S	Gellhorn Short Stem #2, 51mm, 2 in	FVC1615
050157S	Gellhorn Short Stem #3, 57 mm, 2.2 in	FVC1706
050158S	Gellhorn Short Stem #4, 63 mm, 2.5 in	FVC1795
050159S	Gellhorn Short Stem #5, 70 mm, 2.8 in	FVC1676
050160S	Gellhorn Short Stem #6, 76 mm, 3 in	FVC1717
050161S	Gellhorn Short Stem #7, 82mm, 3.2 in	FVC1611
050162S	Gellhorn Short Stem #8, 89mm, 3.5 in	FVC1636

Donut

<u>SKU</u>	<u>Description</u>	<u>SC NPC</u>
050010	Donut #0, 51mm, 2 in	FVC1641
050011	Donut #1, 57mm, 2.2 in	FVC1682
050012	Donut #2, 64mm, 2.5 in	FVC1679
050013	Donut #3, 69mm, 2.7 in	FVC1766
050014	Donut #4, 76mm, 3 in	FVC1718
050015	Donut #5, 82mm, 3.2 in	FVC1597
050016	Donut #6, 88mm, 3.5 in	FVC1666
050199	Donut #7, 95mm, 3.7 in	FVC1610

Shaatz

<u>SKU</u>	<u>Description</u>	<u>SC NPC</u>
050120	Shaatz w/ drain #0 , 38mm, 1.5 in	FVC1730
050121	Shaatz w/ drain #1, 44mm, 1.75 in	FVC1614
050122	Shaatz w/ drain #2, 51mm, 2 in	FVC1659
050123	Shaatz w/ drain #3, 57mm, 2.2 in	FVC1623
050124	Shaatz w/ drain #4, 64mm, 2.5 in	FVC1633
050125	Shaatz w/ drain #5, 70mm, 2.8 in	FVC1760
050126	Shaatz w/ drain #6, 76mm, 3 in	FVC1755
050127	Shaatz w/ drain #7, 82mm, 3.2 in	FVC1729
050128	Shaatz w/ drain #8, 88mm, 3.5 in	FVC1672

Cube

<u>SKU</u>	<u>Description</u>	<u>SC NPC</u>
050066	Cube w/ drain #0, 25mm, 0.98 in	FVC1693
050067	Cube w/ drain #1, 29mm, 1.1 in	FVC1736
050068	Cube w/ drain #2, 32mm, 1.2 in	FVC1697
050069	Cube w/ drain #3, 36mm, 1.4 in	FVC1787
050070	Cube w/ drain #4, 40mm, 1.6 in	FVC1652
050071	Cube w/ drain #5, 45mm, 1.8 in	FVC1808
050072	Cube w/ drain #6, 49mm, 1.9in	FVC1785
050073	Cube w/ drain #7, 55mm, 2.2 in	FVC1635
050074	Cube w/ drain #8, 63mm, 2.5 in	FVC1769
050075	Cube w/ drain #9, 69mm, 2.7 in	FVC1754
050076	Cube w/ drain #10, 74mm, 2.9 in	FVC1707
050055	Cube w/o drain #0, 25mm, 0.98 in	FVC1786
050056	Cube w/o drain #1, 29mm, 1.1 in	FVC1712
050057	Cube w/o drain #2, 32mm, 1.2 in	FVC1774
050058	Cube w/o drain #3, 36mm, 1.4 in	FVC1700
050059	Cube w/o drain #4, 40mm, 1.6 in	FVC1801
050060	Cube w/o drain #5, 45mm, 1.8 in	FVC1713
050061	Cube w/o drain #6, 49mm, 1.9 in	FVC1684
050062	Cube w/o drain #7, 55mm, 2.2 in	FVC1753
050063	Cube w/o drain #8, 63mm, 2.5 in	FVC1616
050064	Cube w/o drain #9, 69mm, 2.7in	FVC1791
050065	Cube w/o drain #10, 74mm, 2.9 in	FVC1680

Cup

<u>SKU</u>	<u>Description</u>	<u>SC NPC</u>
050171	Cup w/ support #0, 48mm, 1.8 in	FVC1678
050172	Cup w/ support #1, 52mm, 2 in	FVC1806
050173	Cup w/ support #2, 57mm, 2.25 in	FVC1756
050174	Cup w/ support #3, 63mm, 2.5 in	FVC1778
050175	Cup w/ support #4, 68mm, 2.7 in	FVC1634
050176	Cup w/ support #5, 73mm, 2.9 in	FVC1719
050177	Cup w/ support #6, 77mm, 3 in	FVC1695
050178	Cup w/ support #7, 82mm, 3.2 in	FVC1607
050163	Cup w/o support #0, 48mm, 1.8 in	FVC1629
050164	Cup w/o support #1, 52mm, 2 in	FVC1677
050165	Cup w/o support #2, 57mm, 2.25 in	FVC1799
050166	Cup w/o support #3, 63mm, 2.5 in	FVC1612
050167	Cup w/o support #4, 68mm, 2.7 in	FVC1657
050168	Cup w/o support #5, 73mm, 2.9 in	FVC1740
050169	Cup w/o support #6, 77mm, 3 in	FVC1622
050170	Cup w/o support #7, 82mm, 3.2 in	FVC1708

Dish

<u>SKU</u>	<u>Description</u>	<u>SC NPC</u>
050137	Dish w/ support #0, 50mm, 2 in	FVC1762
050138	Dish w/ support #1, 55mm, 2.2 in	FVC1810
050139	Dish w/ support #2, 60mm, 2.4 in	FVC1617
050140	Dish w/ support #3, 63mm, 2.5 in	FVC1742
050141	Dish w/ support #4, 68mm, 2.7 in	FVC1703
050142	Dish w/ support #5, 73mm, 2.9 in	FVC1748
050143	Dish w/ support #6, 80mm, 3.1 in	FVC1664
050144	Dish w/ support #7, 85mm, 3.3 in	FVC1609
050129	Dish w/o support #0, 50mm, 2 in	FVC1667
050130	Dish w/o support #1, 55mm, 2.2 in	FVC1654
050131	Dish w/o support #2, 60mm, 2.4 in	FVC1696
050132	Dish w/o support #3, 63mm, 2.5 in	FVC1779
050133	Dish w/o support #4, 68mm, 2.7 in	FVC1775
050134	Dish w/o support #5, 73mm, 2.9 in	FVC1648
050135	Dish w/o support #6, 80mm, 3.1 in	FVC1803
050136	Dish w/o support #7, 85mm, 3.3 in	FVC1764

Oval

SKU	Description	SC NPC
050111	Oval w/ support #1, 51 mm, 2 in	FVC1792
050112	Oval w/ support #2, 57 mm, 2.25 in	FVC1637
050113	Oval w/ support #3, 64 mm, 2.5 in	FVC1710
050114	Oval w/ support #4, 69 mm, 2.7 in	FVC1647
050115	Oval w/ support #5, 76 mm, 3 in	FVC1720
050116	Oval w/ support #6, 83 mm, 3.25 in	FVC1759
050117	Oval w/ support #7, 89 mm, 3.5 in	FVC1665
050118	Oval w/ support #8, 95 mm, 3.75 in	FVC1631
050119	Oval w/ support #9, 102mm, 4 in	FVC1632
050101	Oval w/o support #1, 51mm, 2 in	FVC1626
050102	Oval w/o support #2, 57 mm, 2.25 in	FVC1671
050103	Oval w/o support #3, 64mm, 2.5 in	FVC1658
050104	Oval w/o support #4, 69mm, 2.7 in	FVC1735
050105	Oval w/o support #5, 76 mm, 3 in	FVC1674
050106	Oval w/o support #6, 83 mm, 3.25 in	FVC1765
050107	Oval w/o support #7, 89 mm, 3.5 in	FVC1603
050108	Oval w/o support #8, 95 mm, 3.75 in	FVC1750
050109	Oval w/o support #9, 102 m, 4 in	FVC1687

Marland pessary

SKU	Description	SC NPC
050084	Marland w/ support #2, 57mm, 2.25 in	FVC1781
050085	Marland w/ support #3, 63mm, 2.5 in	FVC1771
050086	Marland w/ support #4, 70mm, 2.75 in	FVC1673
050087	Marland w/ support #5, 76mm, 3 in	FVC1789
050088	Marland w/ support #6, 83mm, 3.25 in	FVC1686
050089	Marland w/ support #7, 89mm, 3.5 in	FVC1773
050090	Marland w/ support #8, 95mm, 3.75 in	FVC1660
050077	Marland w/o support #2, 57mm, 2.25 in	FVC1675
050078	Marland w/o support #3, 63mm, 2.5 in	FVC1699
050079	Marland w/o support #4, 70mm, 2.75 in	FVC1698
050080	Marland w/o support #5, 76mm, 3 in	FVC1797
050081	Marland w/o support #6, 83mm, 3.25 in	FVC1681
050082	Marland w/o support #7, 89mm, 3.5 in	FVC1805
050083	Marland w/o support #8, 95mm, 3.75 in	FVC1749

Hodge pessary

SKU	Description	SC NPC
050045	Hodge w/ support #0, 65mm, 2.6 in	FVC1691
050046	Hodge w/ support #1, 70mm, 2.8 in	FVC1783
050047	Hodge w/ support #2, 75mm, 3 in	FVC1723
050048	Hodge w/ support #3, 80mm, 3.2 in	FVC1655
050049	Hodge w/ support #4, 85mm, 3.4 in	FVC1800
050050	Hodge w/ support #5, 90mm, 3.6 in	FVC1788
050051	Hodge w/ support #6, 95mm, 3.8 in	FVC1690
050052	Hodge w/ support #7, 100mm, 4 in	FVC1751
050053	Hodge w/ support #8, 105mm, 4.2 in	FVC1704
050054	Hodge w/ support #9, 110mm, 4.4 in	FVC1688
050035	Hodge w/o support #0, 65mm, 2.6 in	FVC1722
050036	Hodge w/o support #1, 70mm, 2.8 in	FVC1747
050037	Hodge w/o support #2, 75mm, 3 in	FVC1685
050038	Hodge w/o support #3, 80mm, 3.2 in	FVC1602
050039	Hodge w/o support #4, 85mm, 3.4 in	FVC1731
050040	Hodge w/o support #5, 90mm, 3.6 in	FVC1726
050041	Hodge w/o support #6, 95mm, 3.8 in	FVC1653
050042	Hodge w/o support #7, 100mm, 4 in	FVC1649
050043	Hodge w/o support #8, 105mm, 4.2 in	FVC1752
050044	Hodge w/o support #9, 110mm, 4.4 in	FVC1804

Gehrung pessary

SKU	Description	SC NPC
050145	Gehrung w/ support #0, 50 mm, 1.9 in	FVC1662
050146	Gehrung w/ support #1, 55 mm, 2.2 in	FVC1630
050147	Gehrung w/ support #2, 60 mm, 2.4 in	FVC1782
050148	Gehrung w/ support #3, 65 mm, 2.6 in	FVC1741
050149	Gehrung w/ support #4, 70 mm, 2.8 in	FVC1642
050150	Gehrung w/ support #5, 75 mm, 3 in	FVC1644
050151	Gehrung w/ support #6, 80 mm, 3.2 in	FVC1743
050152	Gehrung w/ support #7, 85 mm, 3.4 in	FVC1724
050153	Gehrung w/ support #8, 90 mm, 3.6 in	FVC1628

TRUSTED PRODUCTS FOR PREMIER WOMEN'S HEALTHCARE

RB Medical Engineering Ltd
Unit 2
Alton Road Industrial Estate
Ross-on-Wye, Herefordshire
HR9 5NS
Tel: +44 (0) 1989 563958
Email: sales@rbmedical.co.uk
Web: www.rbmedical.co.uk

