



Providing the instruments for good health

www.rbmedical.co.uk

Tibbs Cannula

Single Use

SE50899-1	Tibbs Cannula Micro Olive 2.5mm tip (pack of 1) Min 10
SE50900	Tibbs Cannula Olive 3mm tip (pack of 2) Min 10
SE50900-1	Tibbs Cannula Olive 3mm tip (pack of 1) Min 10
SE50901	Tibbs Cannula Kidney/Olive tip (pack of 2) Min 10
SE50901-1	Tibbs Cannula Kidney/Olive tip (pack of 1) Min 10
SE50902	Tibbs Cannula 3mm/5mm Cone tip (pack of 2) Min 10
SE50902-1	Tibbs Cannula 3mm/5mm Cone tip (pack of 1) Min 10
SE50903	Tibbs Cannula 4mm/6mm Cone tip (pack of 2) Min 10
SE50903-1	Tibbs Cannula 4mm/6mm Cone tip (pack of 1) Min 10
SE509031-1	Tibbs Cannula 5mm/8mm Cone tip (pack of 1) Min 10
SE509032-1	Tibbs Cannula 6mm/12mm Cone tip (pack of 1) Min 10





Providing the instruments for good health

www.rbmedical.co.uk

These are used during surgical procedures that require access to veins or arteries, including vascular surgery for irrigation of the vessels to flush out any debris. The primary function of the device is to flush or apply suction to keep the vessels clear of debris during surgery.

Key benefits include:

- Ergonomic tip shapes
- Compatible with all flushing solutions
- Sterile
- Packed individually or as pairs
- High quality
- Easy to use
- Purchase the tip type and size required in the amounts you need (less waste)

Our Tibbs cannula come in a variety of shapes and sizes. If you are unsure of your requirements, please contact your Area Development Manager and they will be able to advise you based on your individual needs.

As a company we have the most exacting standards and comply to the safety and quality standards ISO 13485 and ISO 9001.

Our Tibbs cannula provide a high quality cost effective solution for your surgical needs - if you wish to discuss options available for this product please see our website or contact your Area Development Manager.

- Complies to safety regulations ISO 13485 & ISO 9001
- UK engineered to exacting standards
- Convenient and easy to use