




Providing the instruments for good health

www.rbmedical.co.uk

Easy-part Laparoscopy Instruments

Re-usable

Celebrating 100 years of UK engineering



Key Benefits:

- Convenient and easy to use
- 3 part system - handle, tube & insert
- Large range of inserts
- Extensive scope for combinations
- All parts are compatible & interchangeable

Our easy-part laparoscopy system helps complete diagnosis of a range of abdominal problems. With our extensive range of inserts all needs during surgery can easily be met, whether that be biopsy, grasping or cutting.

The ease and manoeuvrability of our easy-part system can reduce patient trauma and aid patient recovery.

If you wish to discuss options available for this product please see our website or contact your Area Development Manager.

As a company we have the most exacting standards and comply to the safety and quality standards ISO 13485 and ISO 9001 and MDD.

- Complies to safety regulations ISO 13485 & ISO 9001
- UK engineered to exacting standards
- Complies to the Medical Device Directive (MDD)



Handles & Tubes

This is just a sample of the handle options available

www.rbmedical.co.uk



SW10564-10
Plastic Handle
without ratchet



SW10566-10
Plastic Handle
with ratchet



SW10578
Comfort Handle
with ratchet



SW10572
Metal Handle
with ratchet



SW10568
Axial Handle



SW10565
Plastic Handle



SW10554
Insulated Tube
330mm x Ø 10mm

SW10552
Insulated Tube
330mm x Ø 5mm





S5101
5mm Ø, Claw Grasping Forceps 2/3 teeth



S5103
5mm Ø, Double Action Short Forceps with Spoon



S5104
5mm Ø, Fine Dissecting Forceps 1/2 teeth



S5105
5mm Ø, Cobra Forceps 2/4 teeth



S5106
5mm Ø, Spoon Jaw Forceps with Pin



S5107
5mm Ø, Spoon Jaw Forceps without Pin



S5109
5mm Ø, Atraumatic Smooth Jaw Forceps



S5110
5mm Ø, Grooved Forceps with Spoon - Long



S5111

5mm Ø, Grooved Forceps with Spoon - Short



S5112

5mm Ø, Grooved Forceps no Spoon



S5113

5mm Ø, Fine Grooved Forceps 2/3 teeth




S5115

5mm Ø, Serrated Alligator Forceps



S5116

5mm Ø, Fenestrated Blakesley Biopsy Forceps



S5117

5mm Ø, Through Cut Biopsy Punch




S5118

5mm Ø, Universal Double Action Forceps



S5119

5mm Ø, Fine Cross Serration Dissecting Forceps



S5121

5mm Ø, Fine Straight Serration Dissecting Forceps



S5122

5mm Ø, Fenestrated Multi Serration Forceps



S5123

5mm Ø, Mixer Forceps



S5124

5mm Ø, Delphin Forceps



S5125

5mm Ø, Manhes, Wavy Serration Forceps 2/3 teeth



S5126

5mm Ø, Wavy Serration Forceps with Spoon



S5128

5mm Ø, Fine Groove Forceps - Inner Notch



S5129

5mm Ø, Atraumatic Double Action Longitudinal Serrated Forceps



S5130
5mm Ø, Allis Forceps



S5131
5mm Ø, Babcock Forceps



S5132
5mm Ø, Babcock Forceps



S5133
5mm Ø, Mini Fenestrated Forceps



S5134
5mm Ø, Micro Forceps



S5135
5mm Ø, Riddick-Olsen Forceps



S5136
5mm Ø, Riddick-Olsen Forceps



S5137
5mm Ø, De Bakey Forceps



S5138
5mm Ø, De Bakey Forceps



S5139
5mm Ø, Maryland De Bakey Forceps



S5142
5mm Ø, Maryland Mod Forceps



S5143
5mm Ø, Maryland Forceps



S5148
5mm Ø, Johann Fenestrated Forceps



S5149
5mm Ø, Johann Fenestrated Forceps - Long




S5150
5mm Ø, Straight Blunt/Blunt Scissors



S5151
5mm Ø, Curved Metzenbaum Scissors



S5152
5mm Ø, Curved Super Cut Metzenbaum Scissors



S5161
5mm Ø, Serrated Straight Scissors



S5165
5mm Ø, Hooked Scissors



S5166
5mm Ø, Micro Scissors



S5167
5mm Ø, Micro Curved Scissors



S5171
5mm Ø, Left or Right Needle Holders



S5172
5mm Ø, Straight Needle Holders



S5173
5mm Ø, Tungsten Carbide Needle Holders




S5202
10mm Ø, Long Jaw Grasping Forceps



S5205
10mm Ø, Maryland Grasping Forceps



S5207
10mm Ø, Mixer Grasping Forceps



S5208
10mm Ø, Fine Grasping Forceps with Spoon



S5213
10mm Ø, Curved Grasping Forceps - Long



S5214
10mm Ø, Fenestrated Forceps



S5215
10mm Ø, Serrated Alligator Forceps



S5220
10mm Ø, Fenestrated Micro Forceps



S5230
10mm Ø, Square End Serrated Grasping Forceps



S5231
10mm Ø, Pennington Jaw Grasping Forceps



S5232
10mm Ø, Dorsey Jaw Forceps - 30mm



S5233
10mm Ø, Dorsey Jaw Forceps - 40mm




S5234
10mm Ø, Dorsey Jaw Forceps - 50mm




S5235
10mm Ø, Duval Forceps




S5240
10mm Ø, Babcock Forceps



S5243
10mm Ø, Allis De Bakey Forceps



S5244
10mm Ø, Allis Forceps



S5245
10mm Ø, 3 Groove Alligator Forceps



S5246
10mm Ø, Anvil Forceps



S5255
10mm Ø, De Bakey Forceps - 30mm



S5256
10mm Ø, De Bakey Forceps - 40mm




S5257
10mm Ø, De Bakey Forceps - 50mm



S5258
10mm Ø, De Bakey Forceps - 60mm



S5270
10mm Ø, 2/3 Claw Grasping Forceps



S5272
10mm Ø, Biopsy Cup