



Ophthalmic  
Instruments



rbmedical



Re-Usable  
*Catalogue*

[www.rbmedical.co.uk](http://www.rbmedical.co.uk)

# rbmedical

Catalog Number Prefix	Product Group	Page Number
<b>RB10</b>	Speculums & Retractors	04-11
<b>RB20</b>	Knives, Spuds, Blade Breakers, Punches & Handles	12-14
<b>RB30</b>	Lasik/Lasek/Dalk, Curettes, Spatulas, Choppers, Manipulators, Loops, Spoons, Hooks, Fixation Rings, Muscle Hooks, Corneal Markers & Calipers.	16-48
<b>RB40</b>	Scissors	50- 60
<b>RB50</b>	Forceps	62- 100
<b>RB60</b>	Needle Holders	102- 106
<b>RB70</b>	Lacrimal Instruments, Cannulas, Quickswitch I/A System Infusion Handles Bimanual Hand Pieces, Rongeurs, Nasal Speculums, Elevators, Chisels, Mallets & Probes.	107- 113
<b>RB80</b>	Bipolar Forceps.	114
<b>RB90</b>	Sterilization Cases & Trays, Diagnostic Instruments	115-118



# RB10


- Speculums
- Retractors

# RB20

- Knives
- Spuds
- Blade Breakers
- Punches
- Handles

### Barraquer Wire Speculums


**RB10-5000** standard wire 0.8mm, 10mm blades  
**RB10-5009** heavy wire 1mm, 15mm blades  
**RB10-5010** standard wire 0.8mm, 15mm blades



### Barraquer Wire Speculums

Wire Speculum  
**RB10-5355** 10mm blades ( infant )  
**RB10-5356** 15mm blades

Grandon-Barraquer  
**RB10-5011** 15mm blades



### Temporal Approach


Wire Speculum  
 with spring on nasal side  
**RB10-5012** 15mm closed blades  
**RB10-5015** 15mm open blades



### Kratz-Barraquer


Wire Speculum  
 with open blades  
 For unobstructed access to the  
 cataract incision with phaco or I/A tip.

**RB10-5049** heavy wire 1mm, 15 mm blades  
**RB10-5050** standard wire 0.8mm, 15mm blades  
**RB10-5051** standard wire 0.8mm, 13mm blades  
**RB10-18005/A** standard wire 0.8mm, 15mm blades  
**RB10-18005/B** heavy wire 1mm, 15 mm blades



### Barraquer Solid Blade

Wire Speculum  
**RB10-5020** Barraquer, 10mm blades  
**RB10-5030** Barraquer, 15mm blades  
**RB10-5031** Feaster, 20mm blades



## Speculums

### Sauer

Lid Speculum

**RB10-5300** for infants, 11mm blades


**RB10-5302** for premature infants, 2mm blades



### Frankel

Infant Lid Speculum adjustable for precise control of infant eyelid retraction.

**RB10-5310** 10mm blades




### Infant Wire Speculum

**RB10-5320** Blades size 5mm, Open

**RB10-5325** Blades size 5mm, Close

**RB10-5330** Blades size 6mm, Solid



### Alphonso


Infant Lid Speculum finger grips and interlacing blades to facilitate insertion one blade 5mm, one blade 7mm

**RB10-5340**

### Infant Lid Speculum

Lid Speculum closed wire blades gently curved spring

**RB10-5350** 4mm blades



### Cook Eye Speculum

Eye Speculum with locking screw

**RB10-5401** #1, 8mm blades


**RB10-5402** #2, 11mm blades

**RB10-5403** #3, 15mm blades



### Murdoch

Eye Speculum with locking screw  
**RB10-5406** #2, 12mm blades  
**RB10-5407** #3, 15mm blades



### Mellinger

Eye Speculum self retaining  
for unobstructed access to the  
cataract incision with phaco or I/A tip.  
**RB10-5411** 15mm blades



### Williams

Eye Speculum with locking screw  
**RB10-5500** 11mm blades  
**RB10-5510** 15mm blades




### Knapp

Eye Speculum  
**RB10-5511** 13mm blades



### Lancaster Speculum

Eye Speculum  
adjustable mechanism  
with locking nut  
**RB10-5600** 16mm blades



## Speculums

### Lester - Burch Speculum


Eye Speculum with flanges on blades for retaining eye lashes  
adjustable mechanism  
with locking screw  
**RB10-5620**




### Lieberman Speculum

adjustable mechanism


- RB10-5665** Lieberman, K-wire, 5mm blades, infant
- RB10-5670** Lieberman, K-wire, 10mm blades, small
- RB10-5671** Lieberman, K-wire, 15mm blades, standard
- RB10-5675** Lieberman, V-shaped, 15mm blades
- RB10-5668** Tanna-Lieberman, 15mm blades
- RB10-5673** Thorlakson, rounded, 13mm blades
- RB10-5674** Lieberman, triple post, 13mm blades
- RBT10-5671** Titanium 15mm K-wire, blades
- RBT10-5675** Titanium 15mm W-wire, blades
- RBT10-5673** Titanium rounded, 13mm blades




RB10-5670




RB10-5671



RB10-5675



RB10-5668



RB10-5673



RB10-5674

### Lieberman Speculum

Eye Speculum

Adjustable mechanism  
non-reflective finish

**RB10-5666** 5mm blades, Infant

**RB10-5677** 10mm blades, Small

**RB10-5678** 15mm blades, Standard




### Jaffe

Wire Lid Retractor

**RB10-8000** pair, 15mm blades ( large )

**RB10-8001** pair, 13mm blades ( small )



### Castroviejo Speculum

Eye Speculum

adjustable mechanism

**RB10-5700** 5mm blades

**RB10-5701** 10mm blades

**RB10-5702** 14mm blades


**RB10-5703** 16mm blades



### Bangerter Speculum

**RB10-1230** 15mm blades

**RB10-1240** 10mm blades



### Stevenson

Lacrimal Sac Retractor

3x3 blunt prongs, adjustable

**RB10-9500**



### Stevenson

Lacrimal Sac Retractor

solid blades with serrated

edge adjustable

**RB10-9520** 13mm blades





## Speculums

### Weiss



EYE SPECULUM Fenestrated  
**RB10-115950**



### Lieberman Style, aspirating


adjustable mechanism  
 Features aspirating ports to remove  
 excess fluid from the surgical field.  
 Supplied with silicone tubing and  
 luer-lock adaptor.

- RB10-5171** Lieberman, Kratz style, 15mm blades
- RB10-5173** Thorlakson, rounded, 11mm blades
- RB10-5175** Lieberman, V-shaped, 15mm blades
- RB10-5180** Chu, solid, reversible, 15mm blades
- RB10-5179** Replacement silicone tubing, 100cm (39")



### Flieringa

Scleral Fixation Ring  
**RB10-7100** set of eight, 15mm-22mm



### McNeill-Goldman

Scleral Fixation Ring  
 and Blepharostat  
**RB10-7311** Small, 15mm I.D.  
**RB10-7312** Medium, 16mm I.D.  
**RB10-7313** Large, 17mm I.D.




**Schocket**

Scleral Depressor double-end with pocket clip 6mm wide bar and 2.5mm tear drop  
**RB10-7825**



**Josephberg-Besser**

**RB10-7822**



**Schepens**

Orbital Retractor  
**RB10-7880**



**Senn**

Senn Miller Lid Retractor  
**RB10-8100** Sharp, 7x8mm, 16 cm 6¼"  
**RB10-8120** Blunt, 7x8mm, 16 cm 6¼"




**Senn**

Senn Miller Lid Retractor, Double Ended  
**RB10-8150** Sharp  
**RB10-8160** Blunt



**Desmarres**

Lid Retractor thin solid blades  
**RB10-8200**  
 size 0, 11mm wide  
**RB10-8210**  
 size 1, 13mm wide  
**RB10-8220**  
 size 2, 15mm wide  
**RB10-8230**  
 size 3, 17mm wide



## Retractors

### Jaeger

Lid Plate

**RB10-8500** stainless steel, Small

**RB10-8510** stainless steel, Large



### Wright

Fascia Needle

1mm x 6mm oval hole

140mm long

**RB10-8540**




### Knapp

Lacrimal Retractor

4 prongs

**RB10-9000** blunt

**RB10-9001** sharp



### Rollet

Rollet retractor

**RB10-9002** 4-blunt

**RB10-9003** 3-blunt



### Helveston

Tissue Retractor

thin curved blade

**RB10-9010** "Barbie" 7mm wide

**RB10-9012** "Big Barbie" 9mm wide

**RB10-9019** "Great Big Barbie"

11mm wide concave blade



**Tooke**

Corneal Knife  
straight blade  
with curved cutting edge  
**RB20-3650**



**Gill**

Corneal Knife  
curved cutting edge  
**RB20-3700** straight  
**RB20-3710** light curve




**Fukasaku**

"Hockey" Knife  
sharp 3mm front edge  
semi-sharp 6mm  
bottom edge  
**RB20-3735**



**Foreign body spuds**


**RB20-4000** Golf Club, angled, lightly rounded tip  
**RB20-4050** Francis, angled, flag shaped tip  
**RB20-4100** Ellis, lightly curved and rounded tip



RB20-4000




RB20-4050



RB20-4100

**Spud & Gouge**


in reversible handle for  
removing foreign bodies.  
**RB20-4600**



## Punches


### Gass

Corneoscleral Punch  
1.5mm diameter head  
0.75mm deep bite  
rotatable tip for  
convenient positioning  
**RB20-9500**




### Kelly

Descemet's Membrane Punch  
1.0mm diameter head  
0.75mm deep bite  
**RB20-9510**




### Fukasaku

Micro Trabeculectomy Punch  
0.6mm diameter head  
with bullet shaped tip  
0.3mm x 0.6mm deep bite  
**RB20-9512**



### Christensen


Punctum Punch  
1mm wide bite  
for opening stenotic punctums.  
**RB20-9525**



Handles/Blade  
Breakers

**Castroviejo**

Blade breaker and holder  
with sliding lock.  
**RB20-5000**



**Swiss Model**

Blade breaker and holder  
spherical tip.  
**RB20-5150**




**Chuck Handle**

round, serrated  
**RB20-1910** 10cm  
**RB20-1920** 7.5cm



**Bard-Parker Handle**

for disposable slot blades  
**RB20-1950** #3 with gauge  
**RB20-1960** #3



**Barron Handle # 3**

**RB20-1971**



# RB30

- Currettes
- Spatulas/Choppers/
- Manipulators/Loops
- Loops/Spoons/Manipulators
- Hooks/Manipulators
- Fixation Rings
- Muscle Hooks
- Corneal Markers
- Calipers

**Meyerhoefer**

Chalazion Curette  
**RB30-1000** size 0, 1.5mm dia.  
**RB30-1010** size 1, 1.75mm dia.  
**RB30-1020** size 2, 2.25mm dia.  
**RB30-1030** size 3, 2.5mm dia.  
**RB30-1040** size 4, 3.5mm dia.




**Double End Currette**

**RB30-1041** size 1.5mm / 2.5mm dia.  
**RB30-1042** size 2mm / 3mm dia.  
**RB30-1043** size 2.5mm / 3.5mm dia.



**Capsule**

Polishing Curette  
 1.75mm diameter open ring  
**RB30-1100** Nightingale, semi-sharp, posterior edge  
**RB30-1102** Shepherd, semi-sharp, anterior and posterior edges



**Culler**

Iris Spatula  
**RB30-2100** 1mm wide  
**RB30-2110** 2mm wide  
**RB30-2111** nettleship iris repositor



**Wheeler**

Double-End Spatula  
 lightly curved  
**RB30-2120**



**Barraquer**

Iris Spatula  
 round, smooth  
**RB30-2300** 0.5mm diameter  
**RB30-2310** 0.25mm diameter





## Lasik/Lasek/Dalk

### Buratto

LASIK Flap Protector  
For retracting and  
protecting nasally hinged  
corneal flaps.  
**RB30-1770**



### Buratto

LASIK Flap Protector  
blades angled 15° For retracting  
and protecting superiorly  
hinged corneal flaps.  
**RB30-1772**



### Mannis - Buratto

LASIK Flap Protector  
blades angled 45°  
For retracting and protecting  
superiorly hinged corneal flaps.  
**RB30-1773**



### Rowen


Combination LASIK Flap  
Protector and Spatula  
**RB30-1774**



### Dalk

Corneal Dissector  
**RB30-1870**  
12mm long, flattened,  
vaulted spatula with  
beveled tip  
**RB30-1875**


Rosenwasser DALK Trisector  
with triangular shaft and  
spatulated tip



RB30-1870




RB30-1875




**Endothelial Stripper**

1mm wide, thin beveled tip  
**RB30-1830** smooth surface  
**RB30-1835** sosewasser with textured surface



**Sloane**

LASEK Flap Repositor  
 gently curved spatulated tip  
 with beveled smooth edges  
 for repositing and smoothing  
 the epithelial flap after laser  
 ablation  
**RB30-1845**



**Camellin style**

LASEK Trepine &  
 alcohol retaining well  
**RB20-7830** 8mm trephine and 8.5mm well




**Camellin Style**

LASEK Alcohol Well  
 For retaining alcohol solution  
 on the corneal epithelium.  
**RB30-1820** 8.5mm diameter  
**RBT30-1820** 8.5mm diameter titanium  
**RB30-1822** 9.5mm diameter  
**RBT30-1822** 9.5mm diameter titanium



**Anis**


Capsule Scrubber flexible  
 wire textured disc  
**RB30-1111**



## Lasik/Lasek/Dalk


### Moran

LASIK Enhancement  
Spatula Fukasaku spatula  
and a modified paton spatula  
**RB30-2524**



### Corneal Dissector


for intrastromal dissection.  
**RB20-3660** angled, straight blade  
**RB20-3661** angled, curved blade



RB20-3660 RB20-3661

### Fogla Pointed Dissector

Fogla Pointed Dissector  
**RB30-4924**




### Carlson DSEK Smoother

Carlson DSEK Smoother  
**RB30-2040**





### Fogla Trifacet Spatula

Fogla Trifacet Spatula  
**RB30-0479**




**Hersh**

Flap spatula beveled tips  
a combination of a 6mm  
spatula and 10mm spatula  
**RB30-2526**




**Binder**

Flap elevator 1mm wide 12mm  
semi-sharp edges  
**RB30-2528**




**Fukasaku**

LASIK Spatula rounded anterior  
and flat posterior surface with  
beveled tip for lifting and  
dissecting the corneal flap.  
**RB30-2530**




**Guimaraes**

LASIK spatula  
gently curved, bi-convex  
semi-sharp edges  
**RB30-2532**



**Katzen**

LASIK Flap "Unzipper"  
combines a modified  
sinskey hook and a  
bi-convex spatula  
**RB30-2534**



**Maddox**

LASIK Spatula  
combines a flat spatula and  
a cylindrical spatula for  
elevating and repositing the flap  
**RB30-2535**




## Lasik/Lasek/Dalk

### Thorlakson

Flap Lifter combines a modified Sinsky hook and a bi-convex spatula with beveled notches.


**RB30-2536** long handle  
**RB30-2537** short handle



### Lasik Flap Lifter

Combines a 45° angled Sinsky hook and a thin hook with a tip

**RB30-2538**  
**RBT30-2538** Titanium




## Spatulas/Choppers/Manipulators

### Hirschman

Lens spatula flat 0.5mm wide with multiple notches


**RB30-2550**



### Jaffe

Lens spatula flat 0.5mm wide

**RB30-2600**



### Lester

IOL Manipulator  
0.25mm diameter tip ideal for multi-directional manipulation  
**RB30-2680** straight, with guard



IOL Manipulator  
0.25mm diameter tip  
vaulted shank  
with guard  
**RB30-2692**

### Phaco chopper

Chopper 10mm angle to tip  
blunt cutting tip


**RBT30-2180.8**

Chopper 10mm angle to tip  
rounded front, 90° cutting edge.

**RBT30-2180.25**

Delicate 10mm angle to tip  
short tip, 90° cutting edge.

**RBT30-2180.29**



RBT30-2180.8 RBT30-2180.25 RBT30-2180.29



### Mics Phaco Pre-chopper

For bimanual nucleus pre-chopping

**RB30-2323** Alio-Rosen,  
with two defined edges.

**RB30-2324** Alio-Scimitar, curved with  
blunt tip and sharp inferior edge.

### Milder

Phaco Chopper

60° offset inferior wedge  
for horizontal chopping and a  
beveled distal tip for vertical chopping

**RB30-2343** for use in left hand

**RBT30-2343** for use in left hand titanium



RB30-2323 RB30-2324 RB30-2343




### Brierley

Nucleus Splitter

with incision sealing shaft

**RB30-2380**

**RBT30-2380** titanium



### Dodick - Kammann

Nucleus Chopper


for bimanual nucleus division.

**RB30-2382** 1.5mm long tip, str

**RB30-2383** 1.5mm long tip, angled

**RB30-2384** 2mm long tip, straight

**RB30-2385** 2mm long tip, angled



RB30-2382 RB30-2384 RB30-2387

### Escaf

Nucleus Divider

1.25mm long wedge shaped tip  
for vertical chopping.

**RB30-2387**

**RBT30-2387** titanium



## Choppers/Nucleus Manipulators

### Claw Chopper

1.5mm claw shape

**RB30-1446**

### Green nucleus divider

cutting 45 degree to access, 65 degree angle shaft, tip to angle length 10mm

**RB30-1447**

### Ball tip 0.7mm inferior wedge

**RB30-2327**

### Nagahara Phaco Chopper

Rounded, blunt tip 0.9mm inferior wedge

**RB30-2329**

**RBT30-2329 titanium**

### Phaco Chopper

rounded, blunt, olive-tip offset inferior wedge

**RB30-2332** for use in right hand

**RB30-2333** for use in left hand

### Minardi Phaco Chopper

0.7mm long "Quick Chop" tip

**RB30-2370**

**RBT30-2370 titanium**

### Agarwal Karate Phaco Chopper

semi-sharp, wedge shaped inferior edge with a pointed tip.

**RB30-2375**

**RBT30-2375 titanium**

### Chang MicroFinger Chopper

curved finger shaped tip with defined inner edge for horizontal chopping

**RB30-2377**

### Koch Nucleus Chopper

1.6mm long with blunted tip with rounded surfaces designed for safely splitting and dividing a very soft nucleus

**RB30-2386**

**RBT30-2386 titanium**

### Fukasaku wedge shaped inferior edge

**RB30-2392** Snapper Hook

**RB30-2393** Small pupil Snapper Hook with notch to engage the iris

### Rosen Phaco Splitter

wedge-shaped inferior edge blunt tip.

**RB30-2395** inline wedge, standard

**RBT30-2395** inline wedge, standard titanium


**RB30-2396** 60° offset wedge, for use in left hand.



RB30-1446



RB30-1447



RB30-2327




RB30-2329




RB30-2332



RB30-2333




RB30-2370



RB30-2375



RB30-2377



RB30-2386




RB30-2392



RB30-2393



RB30-2395



RB30-2396



RB30-2390


## Nagahara

Nucleus Spatula flat inferior edge

**RB30-2390**

**RBT30-2390 titanium**


**RBM30-2390A** french Pattern



Spatulas/Nucleus Manipulators


**Gimbel**

Nucleus Spatula  
round, smooth  
angled 30° at  
10mm from the tip  
0.25mm diameter  
**RB30-2312**




**Knolle**

Nucleus Spatula  
malleable stainless steel,  
flat 0.5mm wide  
**RB30-2350**




**Koch**

Nucleus Spatula  
gently curved duckbill  
shaped tip  
**RB30-2354** smooth  
**RB30-2355** with notches



**Escaf**


Nucleus Manipulator  
flattened olive-shaped tip angled at 45°  
**RB30-2347**  
**RBT30-2347** titanium



**RB30-2347** **RB30-2357** **RB30-2359**

**Akahoshi**

0.5mm diameter bulbous tip,  
tip length 1.5mm 45° angled  
shaft, tip to angle length 10.0mm  
**RB30-2357**  
**RBT30-2357** titanium



**RB30-2360** **RB30-2361**

**Escaf**

Nucleus Sustainer  
smooth ball tip textured inferior surface  
**RB30-2359**

**Lieberman**

Micro Finger Nucleus Manipulator  
**RB30-2360** for use in right hand  
**RB30-2361** for use in left hand





## Spatulas /Choppers/ Manipulators

### Neuhann

Nucleus Divider  
strongly curved spatulated tip  
**RB30-2388**



### Tennant

Nucleus Rotator  
with 0.6mm diameter  
ball shaped tip  
**RB30-2428**

**RBT30-2428** titanium




### Hunkeler

IOL Manipulator  
with 0.3mm diameter  
hemi-ball tip,  
vaulted shank  
**RB30-2695**

### Mushroom


IOL manipulator angled  
with 0.65mm diameter ball tip  
**RB30-2696**

**RBT30-2696** titanium



### Nucleus sustainer

Nucleus splitter and sustainer,  
curved, blunt tip to be used  
with the prechoppers  
**RBT30-2192.01**




### Scoring reverse hook

"Scoring" reverse hook for DSEK,  
sharp, angled

**RBT30-2183.10**

DSEK stripper,  
"T" shaped tips

**RBT-30-2183.12**



## Spatulas /Choppers/ Manipulators

Double Ended  
Combination of 0.65mm mushroom tip and Nucleus Splitter

**RB30-085**

**RBT30-085 titanium**

Combination of nagahara blunt chopper and sinskey hook

**RB30-086**

**RBT30-086 titanium**

Combination of 0.65mm mushroom and sinskey hook

**RB30-087**

**RBT30-087 titanium**

Combination of nagahara blunt chopper and spatula

**RB30-088**

**RBT30-088 titanium**

Combination of sinskey hook and spatula

**RB30-089**

**RBT30-089 titanium**

Combination of 0.65mm mushroom tip and spatula

**RB30-090**

**RBT30-090 titanium**

Combination of 0.65mm mushroom tip and Y bechert

**RB30-091**

**RBT30-091 titanium**

Combination of sinskey hook and Y bechert

**RB30-092**

**RBT30-092 titanium**

Combination of nagahara blunt chopper and Y bechert

**RB30-093**

**RBT30-093 titanium**

Combination of 0.65mm mushroom tip and sharp chopper

**RB30-094**

**RBT30-094 titanium**

Combination of sinskey hook and sharp chopper

**RB30-095**

**RBT30-095 titanium**

Combination of nagahara blunt chopper and sharp chopper

**RB30-096**

**RBT30-096 titanium**



RBT30-085




RBT30-086



RBT30-087



RBT30-088



RBT30-089



RBT30-090



RB30-091




RB30-092



RB30-093



RB30-094



RB30-095



RB30-096

## Spatulas /Choppers/ Manipulators

Double Ended  
Combination of 0.65mm mushroom tip and Nucleus Splitter

**RB30-085**

**RBT30-085 titanium**

Combination of nagahara blunt chopper and sinskey hook

**RB30-086**

**RBT30-086 titanium**

Combination of 0.65mm mushroom and sinskey hook

**RB30-087**

**RBT30-087 titanium**

Combination of nagahara blunt chopper and spatula

**RB30-088**

**RBT30-088 titanium**

Combination of sinskey hook and spatula

**RB30-089**

**RBT30-089 titanium**

Combination of 0.65mm mushroom tip and spatula

**RB30-090**

**RBT30-090 titanium**

Combination of 0.65mm mushroom tip and Y bechert

**RB30-091**

**RBT30-091 titanium**

Combination of sinskey hook and Y bechert

**RB30-092**

**RBT30-092 titanium**

Combination of nagahara blunt chopper and Y bechert

**RB30-093**

**RBT30-093 titanium**

Combination of 0.65mm mushroom tip and sharp chopper

**RB30-094**

**RBT30-094 titanium**

Combination of sinskey hook and sharp chopper

**RB30-095**

**RBT30-095 titanium**

Combination of nagahara blunt chopper and sharp chopper

**RB30-096**

**RBT30-096 titanium**



RBT30-085




RBT30-086



RBT30-087



RBT30-088



RBT30-089



RBT30-090



RB30-091




RB30-092



RB30-093



RB30-094



RB30-095



RB30-096

**Bechert**

Nucleus Rotator blunt forked tip

**RB30-2400** angled at 7mm

**RBT30-2400** angled at 7mm titanium

**RB30-2410** angled at 9mm

**RBT30-2410** angled at 9mm titanium

**Bechert**

Lens Pusher vertical Y-shaped tip angled.

**RB30-2727**

**RBT30-2727** titanium

**Drysdale**

Nucleus Manipulator paddle shaped tip

**RB30-2418** rounded edges

**RBT30-2418** rounded edges titanium

**RB30-2419** defined anterior edge for capsule polishing

**Pisacano**

Nucleus Rotator flat disc shaped tip

**RB30-2420** horizontal

**RBT30-2420** horizontal titanium

**RB30-2421** vertical

**RBT30-2421** vertical titanium

**Siepser**

Endocapsular Controller trumpet shaped tip

**RB30-2425**

**Anis**

Nucleus Manipulator with 0.5mm diameter ball tip

**RB30-2426** straight

**RB30-2427** angled

**Grayson**

Nucleus Manipulator flattened, smoothly polished 1.1mm long with "tear-drop" shaped tip.

**RB30-2417**

**Lester**

IOL Manipulator 0.25mm diameter tip Ideal for multi-directional manipulation.



**RB30-2690**

**RBT30-2690** titanium

IOL Manipulator 0.25mm diameter tip vaulted shank For engaging the inferior haptic without touching the dome of the IOL.

**RB30-2691**

**RBT30-2691** titanium



## Hooks/Manipulators

### Reverse Sinskey Hook

0.2mm diameter blunt tip for scoring the recipient endothelium and positioning the donor lamella.

**RB30-5002**

### Sinskey

Lens Manipulating Hook round handle  
0.25mm diameter blunt tip

**RB30-5192** straight, with guard

**RB30-5193** angled

**RBT30-5193** angled titanium

### Sinskey II

Lens Manipulating Hook  
0.2mm diameter blunt tip round handle

**RB30-5232** straight, with guard

**RB30-5233** angled

**RBT30-5233** angled titanium

### Lewicky

Lens Manipulating Hook  
0.15mm dia blunt tip ang vaulted shank

**RB30-5250**

**RBT30-5250** titanium

### Hirschman

Iris Hook blunt, angled

**RB30-5400**

### Bonn

Micro Iris Hook blunt, flat handle

**RB30-5060**

### Shepard

Iris Hook blunt

**RB30-5110**

### Lens Manipulating Hook

0.15mm diameter blunt tip

**RB30-5160** Jaffe, straight

**RB30-5170** Clayman-Jaffe, angled

### Sinskey

Lens Manipulating Hook  
0.25mm diameter blunt tip straight

**RB30-5180**

### Harris - Sinskey

Lens Manipulating Hook  
0.25mm diameter blunt tip angled


**RB30-5190**

### Sinskey II

Lens Manipulating Hook  
0.2mm diameter blunt tip

**RB30-5220** straight


**RB30-5230** angled




RB30-5002 RB30-5192 RB30-5193 RB30-5232



RB30-5233 RB30-5250 RB30-5400



RB30-5060 RB30-5110 RB30-5160 RB30-5170




RB30-5180 RB30-5190 RB30-5220 RB30-5230



**Lippman**

Nucleus Spatula  
0.5mm diameter round  
shaft with blunt forked tip  
**RB30-2430**



**Girard**

Synechiaie Spatula  
0.5mm wide with  
semi-sharp tip  
**RB30-2500**



**Obstbaum**

Synechiaie Spatula  
0.25mm wide with  
semi-sharp tip  
**RB30-2520**



**Benedetto**

Spatula 10mm long,  
smoothly polished  
flattened spatula  
with beveled notches  
**RB30-2522**



**Elsching**

Cyclodialysis Spatula  
flat, 1mm wide.  
**RB30-3000**



## Spatulas/ Manipulators

### Castroviejo

Synechiaie Spatula  
0.5mm wide with  
semi-sharp tips.  
**RB30-3150**



### Castroviejo

Cyclodialysis Spatula  
**RB30-3200** 0.25mm wide  
**RB30-3210** 0.5mm wide


### Kimura

Platinum Spatula  
**RB30-3500**



### Nucleus Dividing and Separating


**RB30-52001** size: 0.6×1.2mm tip oval  
**RB30-52003** size: 0.4×0.65mm tip oval  
**RB30-52004** size: 0.4×0.65mm tip oval



### Straight Ring Manipulator

Straight shaft  
0.5mm inside diameter, 0.8mm outside diameter, rounded tip  
Round handle, length 121mm  
**RB30-6462**

0.5mm inside diameter, 0.8mm outside diameter, rounded tip  
Angled shaft, tip to angle length 10mm  
Round handle, length 118mm  
**RB30-6462-1**



0.8mm diameter, rounded tip  
Angled shaft, tip to angle length 10mm  
Round handle, length 118mm  
**RB30-6462-2**



**Keener**

Lens Divider  
with three preformed snare  
wire loops  
**RB30-3900** complete  
**RB30-3910** wire loops



**Lens Loop**

**RB30-4000** Lewis, small  
**RB30-4010** Lewis, Large  
**RB30-4030** New Orleans  
**RB30-4100** Wilder, serrated



**Kirby**

Lens Loop and  
Muscle Hook  
**RB30-4105**



**Morrison**

Lens Loop and  
Probe  
**RB30-4120**



**Lens Loop**

serrated, 3mm wide  
**RB30-4140** Kansas, light curve  
**RB30-4141** Keener-Art, strong curve





**Alfonso**


Cutting Platform  
Used as a cutting  
block for nucleus  
triselector.

**RB30-4156** 2mm wide  
**RB30-4157** 3mm wide




**Giannetti**

Nucleus Capture Spoon  
**RB30-4160**



**Culler**

Lens Spoon  
sterling silver  
**RB30-4200**




**Paton**

Spatula & Spoon  
For corneal  
button transfer  
**RB30-4255**



**Bunge**

Evisceration Spoon  
**RB30-4400** small  
**RB30-4410** large



**Wells**

Enucleation  
Spoon  
**RB30-4450**




**Graether**

Collar Button  
micro iris retractor  
**RB30-4900** straight, with guard  
**RB30-4910** angled



**Keuch**

Small Pupil Dilator  
with guard for enlarging  
a small pupil through the  
sideport and primary incision.  
**RB30-4950**




**ICL Manipulator**

**RB30-5011**



**Bonn**


Micro Iris Hook  
extra delicate blunt,  
round handle with guard.  
**RB30-5050**



## Hooks/Manipulators


### Kuglen Iris Hook & Lens Manipulator

- RB30-5500** Push-pull model, flat handle, straight
- RB30-5520** Push-pull model, flat handle, angled
- RB30-5510** Clover leaf model, flat handle, straight
- RB30-5522** Push-pull model, round handle, straight w/guard
- RB30-5523** Push-pull model, round handle, angled





### Lens Manipulator

- straight, with guard
- RB30-5526** Rentsch, "Boat Hook"
- RB30-5528** Osher, Y-hook




### Lens Manipulating Hook

- 0.15mm diameter tip
- RB30-5530** Fenzl, round handle straight with guard
- RB30-5531** Fenzl, round handle, angled
- RB30-5532** Maltzman-Fenzl, round handle, straight with guard
- RB30-5533** Maltzman-Fenzl, round handle, angled
- RB30-5535** Jaffe-Maltzman, flat handle, angled



### Deltz

- 0.15mm diameter tip
- RB30-5566**



### Alio - Kelman

- IOL Hook 0.2mm diameter tip
- RB30-2697**




**Fixation Hook**

double, sharp  
**RB30-6000**  
**RB30-6010**




**Twist Hook**

for scleral fixation  
**RB30-6100** right  
**RB30-6110** left




**Silverman**

Fixation Pick  
 blunt cone shaped tips  
 spaced 3mm apart  
 for atraumatic scleral fixation  
**RB30-6130**



**Globe Ring**

Globe Fixation Ring  
 with concentric grooves  
**RB30-6153** for refractive procedures  
**RB30-6154** for cataract procedures



**Lu - Mendez**

LRI Guide and Fixation Ring  
 with delicate blunt teeth marks  
 10° increments 11mm I.D.  
**RB30-6158**  
**RBT30-6158** titanium



## Fixation Rings


### Gimbel

Stabilization Ring  
with swivel handle  
13mm diameter with blunt teeth  
on both sides  
For stabilizing either the right or  
left eye during LASIK.  
**RB30-6160**



### Fine - Thornton

Swivel Fixation Ring  
blunt teeth features marks on anterior  
surface in 15° increments as a guide  
for placement of limbal relaxing incisions.  
**RBT30-6161** 13mm diameter titanium  
**RB30-6161** 13mm diameter  
**RB30-6163** 16mm diameter



### Thornton

Globe Fixation Ring  
multiple blunt teeth  
16mm diameter ring  
**RB30-6165**  
**RBT30-6165** titanium



### Hofmann - Thornton

Globe Fixation Ring  
with swivel handle  
16mm diameter ring  
multiple blunt teeth  
**RB30-6168**



**Stevens**

Tenotomy Hook  
blunt, small  
**RB30-6600**




**Helveston**

"Teaser" Hook  
extra delicate, blunt  
**RB30-6610**



**Graefe**

Muscle Hook  
**RB30-6710** size 1  
**RB30-6720** size 2  
**RB30-6730** size 3




**Gass**

Retinal Detachment  
Hook with oval hole  
**RB30-6740**



**Helveston**

Muscle Hook  
extra delicate,  
with blunt spatulated tip  
**RB30-6760** large  
**RB30-6761** small



## Muscle Hooks


### Helveston

Muscle Hook  
extra delicate, with blunt  
spatulated tip  
**RB30-6778** 8mm  
**RB30-6780** 10mm  
**RB30-6782** 12mm



### Jameson

Muscle Hook  
**RB30-6800** small  
**RB30-6810** large



### Du Plessis

Muscle Hook  
double arm  
**RB30-6814**



### Guyton

Small Incision  
Muscle Hook  
**RB30-6820**



### Green

Muscle Hook  
**RB30-6850**



### Moorfields


Moorfields  
Strabismus Hook  
**RB30-6851**



**Mendez**

Degree Gauge 12mm I.D. marked in 10° increments for accurate axial alignment on the cornea

**RBT30-7900 titanium**  
**RB30-7900**



**Bores**

Axis Marker


**RBT30-7910 titanium**  
**RB30-7910**



**Henderson**

Degree Gauge 11.7mm I.D. Beveled face with markings in 5° increments 4 notches on internal edge to align with cardinal positions


**RBT30-7904 titanium**



**Henderson**

Toric IOL Marker two 2.5mm radial marks extending to 11mm (used with ST30-7904)


**RBT30-7912 titanium**



**Henderson**

Henderson Alignment Marker with atraumatic disc-shaped marking surfaces at 3, 6 & 9 o'clock

**RBT30-7908 titanium**





## Corneal Markers

### Optic Zone Markers

#### Bores

- RB30-8130 3.0mm
- RB30-8131 3.25mm
- RB30-8132 3.5mm
- RB30-8134 4.0mm
- RB30-8135 4.25mm
- RB30-8136 4.5mm
- RB30-8137 4.75mm
- RB30-8138 5.0mm
- RB30-8142 6.0mm
- RB30-8146 7.0mm
- RB30-8150 8.0mm
- RB30-8154 9.0mm

#### Hoffer, with cross hairs

- RB30-8330 3.0mm
- RB30-8340 3.25mm
- RB30-8360 3.75mm
- RB30-8370 4.0mm
- RB30-8380 4.25mm
- RB30-8390 4.5mm
- RB30-8400 4.75mm
- RB30-8410 5.0mm
- RB30-8412 5.5mm
- RB30-8414 6.0mm
- RB30-8418 7.0mm
- RB30-8422 8.0mm
- RB30-8426 9.0mm



### Anis

#### Radial Marker

- 4mm diameter center opening
- 12.5mm O.D.
- RB30-8852 8 blades
- RB30-8854 12 blades
- RB30-8855 16 blades



### Anis

#### Suture Placement Marker

- 8 radial blades
- 8mm diameter center opening for corneal grafts.
- RB30-8862



### Fechtner

#### Trabeculectomy Marker

RB30-8975



MARKING PATTERN



**McDonald**

Optic Zone Marker  
double end, with cross hairs  
**RB30-8117**  
6.5 and 7.0mm diameters



**Osher Neumann**

Osher-Neumann  
Corneal Marker  
**RB30-8000** 4 Line  
**RB30-8001** 8 Line  
**RB30-8002** 12 Line



**PRK Scraper**

combination of  
PRK Scraper and Marker  
6.5mm marker  
total length 112mm  
**RB30-8003**



**Double Ended Marker**

5mm diameter on one side  
5.5mm opposite  
**RB30-8004**



**Maloney**

Intra-Operative  
Keratometer titanium  
for qualitative measurement  
of astigmatism.  
**RBT30-7888** titanium  
**RB30-7888**



## Calipers


### Castroviejo

Caliper  
20mm spread  
**RB30-9000** straight  
**RB30-9002** angled  
**RB30-9010** curved




### Helveston

Scleral Ruler  
15mm long, with notches  
in 5mm increments  
single marking tooth at tip for  
measuring and marking the  
distance from the limbus to  
the location of the muscle  
reinsertion.  
**RB30-9030**



### Jameson

Caliper  
chrome plated brass  
**RB30-9150**



### Ruler

Stainless Steel  
inches and mm  
**RB30-9200**



### Braunstein

Caliper  
3.5/4.0mm  
Stainless steel  
**RB30-9250**  
**RB30-9251**



# PRECISE, CONTROLLED DELIVERY, TIME AFTER TIME

## Y-TIP PLUNGER ROD

- Rounded, Capsular-friendly dimensions assist lens manipulation post-implantation

Distinctive blue color designed for easy visibility in the cartridge and reduced microscope reflections  
**RBT50-7791**



## CARTRIDGE LOCK

Quick and easy snap-in design, for secure cartridge loading, assures correct position of cartridge every time.

## SCREW-STYLE

Consistent insertion forceps result in predictable lens delivery from start to finish.

High quality titanium material and durable construction for safe, reliable implantation




rbmedical

# RB40

● Scissors

### Vitreotomy Scissors


“Squeeze Handle” Scissors  
these “squeeze handle” instruments feature titanium handles to reduce overall weight, and front ends of hardened stainless steel for secure cutting



#### Scissors

horizontal, straight


- RB40-5300** 20G 0.8 S.S
- RBT40-5300** 20G 0.8 titanium
- RB40-5301** 23G 0.6 S.S
- RBT40-5301** 23G 0.6 titanium



#### Scissors

horizontal, curved


- RB40-5310** 20G 0.8 S.S
- RBT40-5310** 20G 0.8 titanium
- RB40-5423** 23G 0.6 S.S
- RBT40-5423** 23G 0.6 titanium
- RB40-5425** 25G 0.50 S.S
- RBT40-5425** 25G 0.50 titanium



#### Scissors

horizontal, angled 45°

- RB40-5320** 20G 0.8 S.S
- RBT40-5320** 20G 0.8 titanium
- RB40-5321** 23G 0.6 S.S
- RBT40-5321** 23G 0.6 titanium



#### Scissors

vertical opening straight


- RB40-5200** 20G 0.8 S.S
- RBT40-5200** 20G 0.8 titanium
- RB40-5201** 23G 0.6 S.S
- RBT40-5201** 23G 0.6 titanium



#### Scissors

vertical opening angled 60°

- RB40-5240** 20G 0.8 S.S
- RBT40-5240** 20G 0.8 titanium



#### Scissors

vertical opening angled 80°


- RB40-5250** 20G 0.8 S.S
- RBT40-5250** 20G 0.8 titanium

### Alio


MICS scissors

23 gauge curved shaft with horizontal micro blades for cutting delicate membranes

- RBT40-5351** titanium



Horizontal blade



Vertical blade

### YAE

Multi - Purpose 23G vertical blade

- RBT40-5352** 3mm blade 24mm long shaft 112mm
- RBT40-5353** 2.5mm blade, 15mm long curved shaft 100mm



## Scissors


### Castroviejo

Corneoscleral Scissors  
small blades  
**RB40-1000** right  
**RB40-1010** left



### Castroviejo

Corneoscleral Scissors  
medium blades  
**RB40-1020** right  
**RB40-1030** left



### Troutman - Castroviejo

Corneoscleral Scissors  
with stop, small blades  
**RB40-1200** right  
**RB40-1210** left



### Troutman - Castroviejo

Corneoscleral Scissors  
with stop, medium blades  
**RB40-1220** right  
**RB40-1230** left



### Castroviejo

Corneal Transplant Scissors  
with stop  
strong curve, small blades  
**RB40-2000** right  
**RB40-2010** left




## Retinal Scissors

### Rich

Corneal Transplant Scissors extremely delicate, strongly curved blades for trimming the recipient corneal bed and donor button. Can also be used for open sky anterior capsulotomy.


**RB40-2014** right  
**RB40-2015** left



### Terry

Micro Corneal Scissors round handle for meticulous control.

**RB40-2016** right  
**RB40-2017** left



### Katzin

Corneal Transplant Scissors strong curve, medium blades


**RB40-2020** right  
**RB40-2030** left



### Castroviejo

Universal Corneal Scissors curved, blunt tips

**RB40-2200** small blades  
**RB40-2220** medium blades  
**RB40-2240** large blades  
**RB40-2260** strongly curved large blades



RB40-2200

RB40-2220

RB40-2240

RB40-2260



### Castroviejo

Keratoplasty Scissors angled to side blunt tips, medium blades

**RB40-2850**






## Scissors

### Corneal scissors

curved, 10mm long blades,  
**RB40-1464-1** pointed tips curved  
**RB40-1464-2** blunt tips




RB40-1464-1





### Conjunctival scissors

curved, 10mm long blades,  
**RB40-1540-1** pointed tips curved  
**RB40-1540-2** reverse blades  
**RB40-1540-2** blunt tips




### Conjunctival scissors

curved, 11mm long blades,  
 blunt tips  
**RB40-1543**




### Conjunctival scissors

curved, 8mm long blades,  
 blunt tips  
**RB40-1546**



### IOL removal scissors

Cutting scissors for removal of  
 IOL straight, blunt tips, 8mm long  
 with serrated lower blade  
**RB40-1690**




### IOL removal scissors

Cutting scissors for removal of  
 IOL angled, blunt tips, 8mm long  
 with serrated lower blade  
**RB40-1691**




**McPherson - Westcott**

Conjunctival Scissors curved  
blunt tips, small blades.  
**RB40-3000**



**Westcott**

Tenotomy Scissors  
blunt tips, medium blades  
curved.  
**RB40-3004**




**Westcott**

Tenotomy Scissors  
blunt tips, standard blades curved  
**RB40-3100** right  
**RB40-3110** left




**Tibolt**

Punctal Scissors  
one conically shaped blade  
with micro grooves and  
textured surface for simple  
"Two-snip Punctoplasty"  
**RB40-3130**



**Mini Westcott**

Tenotomy Scissors  
4mm long curved blades  
blunt tips.  
**RB40-3200**



## Scissors

### McPherson - Westcott

Stitch Scissors  
very sharp pointed  
tips small blades curved.  
**RB40-4000**



### Westcott

Stitch Scissors  
very sharp pointed tips  
standard blades curved  
**RB40-4100**



### Jaffe

Stitch Scissors  
very sharp pointed tips  
medium blades curved  
**RB40-4150**



### Iris

Scissors  
sharp tips, small blades  
**RB40-4400** straight  
**RB40-4410** curved



### Fine

Stitch Scissors  
sharp tips, medium blades  
**RB40-4500** straight  
**RB40-4510** curved




**Uthoff**

Capsulotomy scissors,  
convex curved,  
7 mm long, thin blades,  
pointed tips, right  
**RB40-4240**





**Fogla Scissors Right**

The lower blade is thinner and  
1mm longer than the upper blade  
**RB40-4241**




**Multipurpose scissors Noyes**

Multipurpose scissors Noyes,  
straight blades, blunt tips  
**RB40-4242**



**Corneal scissors**

Corneal scissors, strong, angled,  
8 mm long blades, blunt tips  
**RB40-4243**



## Scissors


### Barraquer

Iris Scissors  
blunt tips, 7mm long blades.  
**RB40-4900**



### Capsulotomy Scissors Vannas

5mm long blades, sharp tips  
**RB40-5000** straight  
**RB40-5010** curved  
**RB40-5020** angled forward  
**RB40-5030** angled to side



### Gills - Vannas

7mm long blades, sharp tips  
**RB40-5050** straight  
**RB40-5060** curved  
**RB40-5070** angled forward




### Mcperson - Vannas

8mm long blades, sharp tips.  
**RB40-5100** curved

### Capsulotomy Scissors Vannas

7mm extra thin long blades, sharp tips  
**RB40-5110** straight  
**RB40-5111** curved  
**RB40-5112** angled forward  
10mm extra thin long blades, sharp tips  
**RB40-5120** straight  
**RB40-5121** curved  
**RB40-5122** angled forward  
**RB40-5123** curved, angled, forward  
10mm long angled blades, curved tips  
**RB40-5125** sharp points  
**RB40-5126** blunt points



### Eye Scissors

pointed tips small  
**RB40-7200** straight  
**RB40-7210** curved



Scissors

Eye Scissors

pointed tips  
standard

**RB40-7400** straight

**RB40-7410** curved

**RBTC40-7400** Tungsten Carbide Str

**RBTC40-7410** Tungsten Carbide Cvd

**RBSC40-7400** Super Cut Str

**RBSC40-7410** Super Cut Cvd



Eye Scissors

large rings  
pointed tips  
standard

**RB40-7440** straight


**RB40-7450** curved



Stitch Scissors

ribbon type  
needle points curved

**RB40-8000**



Stitch Scissors

needle points  
light curve

**RB40-8100**



## Scissors

### Stevens

Tenotomy Scissors

blunt tips

standard

**RB40-8500** straight

**RB40-8510** curved

**RBTC40-8500** Tungsten Carbide Str

**RBTC40-8510** Tungsten Carbide cvd

**RBSC40-8500** Super Cut str

**RBSC40-8510** Super Cut cvd




### Stevens

Tenotomy Scissors

ribbon type blunt tips

**RB40-8520** straight

**RB40-8530** curved



### Knapp

Strabismus Scissors

blunt tips

**RB40-9000** Straight


**RB40-9010** curved

**RBTC40-9000** Tungsten Carbide Str

**RBTC40-9010** Tungsten Carbide Cvd

**RBSC40-9000** Super Cut str

**RBSC40-9010** Super Cut cvd




### Strabismus Scissors

ribbon type

blunt tips

**RB40-9020** straight

**RB40-9030** curved



### Metzenbaum

Dissecting Scissors  
gently curved blades  
blunt tips  
5 3/4" (14.5cm)  
long

**RB40-9040**

**RBTC40-9040** Tungsten Carbide

**RBSC40-9040** Super Cut



### Utility Scissors

one pointed and one  
blunt tip  
5 3/4" (14.5cm) long  
straight

**RB40-9060**

**RBTC40-9060** Tungsten Carbide

**RBSC40-9060** Super Cut




### Enucleation Scissors

blunt tips

**RB40-9910** medium curve

**RB40-9920** strong curve



### Black Coated Super Cut Scissors

black coated super cut Iris scissors

**RBSCB40-7400** Str

**RBSCB40-7410** Cvd





rbmedical

# RB50

● Forceps

### ICL forceps

ICL forceps for loading, angled tips  
**RB50-75666** S.S  
**RBT50-75666** titanium




### Foldable Lens Cutter Squeeze

Titanium Handle 112mm  
**RBT-4-21731**



### IOL Lens Cutter

IOL Lens cutter/Holder with one serrated blade stainless steel  
**RBT50-0787**



### Foldable Lens Removing Forceps

Titanium Handle 110mm  
**RBT-4-21741**



# Retinal Forceps

## Vitreo - Retinal Forceps

Stainless Steel Autoclavable




smooth jaws  
20 gauge, straight  
**RB50-7499** 20G/0.9



Stainless Steel Autoclavable



Subretinal Forceps  
Angled 135°,  
4.5mm from the tip  
**RB50-7500** 20G/0.9



Serrated jaws  
straight  
**RB50-7501** 20G/0.9  
**RB50-7502** 23G/0.6



ILM Forceps  
2.2mm long with  
shaft size 27.5mm  
**RB50-7503** 20G/0.9  
**RB50-7504** 23G/0.6



End gripping forceps  
straight  
**RB50-7505** 20G/0.9  
**RB50-7506** 23G/0.6



Serrated ILM  
Forceps, 1.25mm  
Long – with fine pick.  
**RB50-7507** 20G/0.9  
**RB50-7508** 23G/0.6



Spaide Pic forceps,  
**RB50-7510** 20 G/0.9  
**RB50-7511** 23 G/0.6



Asymmetrical  
End-gripping Forceps,  
**RB50-7512** 20 G/0.9  
**RB50-7513** 23 G/0.6



Eckardt End-gripping  
Forceps,  
**RB50-7514** 20 G/0.9  
**RB50-7515** 23 G/0.6



ILM forceps,  
**RB50-7516** 20 G/0.9  
**RB50-7517** 23 G/0.6



Foreign Body Forceps  
Basket style with three  
gripping arms.  
**RB50-7518** 20 G/0.9

Vitreo - Retinal Forceps



Titanium



Subretinal Forceps  
Angled 135°,  
4.5mm from the tip

**RBT50-7500 20G/0.9**




Serrated jaws  
straight

**RBT50-7501 20G/0.9**  
**RBT50-7502 23G/0.6**



ILM Forceps  
2.2mm long with  
shaft size 27.5mm

**RBT50-7503 20G/0.9**  
**RBT50-7504 23G/0.6**



End gripping forceps  
straight

**RBT50-7505 20G/0.9**  
**RBT50-7506 23G/0.6**



Serrated ILM  
Forceps, 1.25mm  
Long – with fine pick.

**RBT50-7507 20G/0.9**  
**RBT50-7508 23G/0.6**




Spaide Pic forceps,

**RBT50-7510 20 G/0.9**  
**RBT50-7511 23 G/0.6**



Asymmetrical  
End-gripping Forceps,

**RBT50-7512 20 G/0.9**  
**RBT50-7513 23 G/0.6**




Eckardt End-gripping  
Forceps,

**RBT50-7514 20 G/0.9**  
**RBT50-7515 23 G/0.6**



ILM forceps,

**RBT50-7516 20 G/0.9**  
**RBT50-7517 23 G/0.6**



Foreign Body Forceps  
Basket style with three  
gripping arms.

**RBT50-7518 20 G/0.9**

## Titanium Retinal Forceps

### Micro Capsulorhexis Forceps

Capsulorhexis forceps 4mm long jaws with angled, grasping platform. Curved shaft 24mm long.

**RB50-7651** 21G/0.8 S.S

**RBT50-7651** 21G/0.8 titanium

**RB50-7652** 23G/0.6 S.S


**RBT50-7652** 23G/0.6 titanium



Ikeda Capsulorhexis forceps  
0.6mm (23G) Delicate 1.2mm long jaws  
with angled, grasping platform.  
Curved shaft 24mm

**RB50-7653** S.S

**RBT50-7653** titanium



Ikeda angled micro Capsulorhexis  
0.9mm long jaw with  
90° angled,  
grasping platform.  
Strongly curved shaft 24mm

**RB50-7654** S.S

**RBT50-7654** titanium

### 23 G, MICS/MCCS Coaxial Capsulorhexis Forceps

MICS/MCCS Capsulorhexis 23G  
1.3 mm jaws, for 1 mm incisions

**RBT50-600001** titanium



### MICS capsulorhexis forceps

**RBT50-7655**

23G Ikeda capsulorhexis forceps 0.65 mm  
delicate 1.2 mm long jaws with angled,  
grasping platform, curved shaft

**RBT50-7656**

23G komaki micro capsulorhexis forceps  
universal capsulorhexis forceps 0.7 mm  
long jaws with 90° angled grasping platform




### Botvin

Iris Forceps  
For peripheral iridectomies.  
**RB50-1040** serrated tips  
**RB50-1050** 1x2 teeth




### Colibri Style

very delicate, 1x2 teeth, 0.12mm  
**RB50-1500** Colibri  
**RB50-1520** Harms-Colibri, with tying platform  
**RBT50-1520** Harms-Colibri, titanium, with tying platform  
**RB50-1522** RByson-Colibri, short tips, extremely thin shanks, with tying platform Standard, 1x2 teeth with tying platform  
**RB50-1550** Colibri, 0.4mm Corneal Suturing Forceps delicate, 1mm spread, 1x2 teeth, 0.12mm  
**RB50-1560** Polack Corneal Forceps very delicate, 1x1 teeth, 0.1mm  
**RB50-1700** Pierse-Colibri




### Gills - Colibri

Tissue Forceps  
very fine pointed tips  
with tying platform  
**RB50-1450**



### Colibri Forceps

Colibri corneal forceps,  
delicate, 1x2 straight,  
with tying platform  
**RB50-118950-1** 0.12 mm  
**RB50-118950-2** 0.40 mm




### Harms

**RB50-118950-3** 0.12 mm teeth,  
with tying 1x2 oblique

### Hoskin

Tissue Forceps, straight  
**RB50-1718** #18, fine grooved tips  
**RB50-1719** #19, micro grooved tips  
**RB50-1728** #28, very delicate, micro grooved tips



## Forceps


### Hoskin

Tissue Forceps,  
Colibri style


**RB50-1722** #22, micro grooved tips

**RB50-1730** #30, very delicate micro  
grooved tips

**RB50-1734** #34, pointed tips



RB 50-1722



RB 50-1730




RB 50-1734



### Pierse

Corneal Forceps  
very delicate 1x1 teeth,  
0.1mm straight


**RB50-1800**



### Fechtner

Conjunctiva Forceps  
very delicate ring shaped  
tips with tying platform

**RB50-1820**




### Castroviejo - Colibri

**RB50-1900** 1x2 teeth, 0.12mm

**RB50-1902** Harms 1x2 teeth,  
0.12mm with tying platform



**RB50-1920** Hyde 2x2 teeth,  
0.12mm with tying platform



### Hoskin

Tissue Forceps, curved

**RB50-1716** #16, fine grooved tips




**Bonn**

French pattern Bonn forceps 0.1 mm with tying platform  
**RB50-7850A** 11cm long handle  
**RB50-7835** 7.5cm short handle  
**RBT50-7835** titanium




**Bonn**

Forceps  
 very delicate 1x2 teeth, 0.12mm  
**RB50-2000**  
**RB50-2010** with tying platform




**Bonn**

Forceps  
 1x2 teeth, 0.12mm  
 with tying platform  
**RB50-2020** S.S  
**RBT50-2020** titanium




**Thorpe**

Corneal Forceps  
 2x3 teeth, 0.12mm  
 with tying platform  
**RB50-2100**



**Colibri Corneal Forceps**

Colibri Corneal Forceps  
 with tying platform titanium  
**RBT50-500095** 0.10 1x2  
**RBT50-500095-1** 0.30 1x2






## Forceps


### Bonn Corneal Forceps

Bonn Corneal Forceps  
delicate, 1x2 straight 0.12mm teeth  
with tying platform  
**RB50-500093**




### Corneal Forceps

Corneal Forceps  
delicate, 1x2 straight 0.12mm teeth  
with tying platform, ergonomic handle  
**RB50-500094**




### Corneal Forceps

Corneal Forceps 1x2 delicate w/tying  
platform titanium ergonomic handle  
**RBT50-500096** 0.10  
**RBT50-500096-1** 0.30  
**RBT50-500096-2** 0.10 oblique  
**RBT50-500096-3** 0.30 oblique



### Corneal Forceps

Corneal Forceps 90mm titanium  
**RBT50-500097** 0.12 1x2  
**RBT50-500097-1** 0.3 1x2  
**RBT50-500097-2** 0.12 1x2 oblique  
**RBT50-500098** 0.3 notched  
Plain tying forceps  
**RBT50-500098-1** 5mm titanium




### Kelmann McPherson

Kelmann Forceps 90mm titanium  
**RBT50-600009** tying platform Ang 9mm  
**RBT50-600009-1** tying platform Ang 7mm  
**RBT50-600009-2** tying platform Ang 11mm




### Iris forceps

Iris forceps, 1x2 straight,  
0.20 mm teeth, with tying platform  
**RB50-700006**




### Corneal forceps

Corneal forceps, delicate, 1x2 straight  
0.12 mm teeth, with tying platform  
**RB50-13160**



### Conjunctival and Fixation Forceps


Conjunctival and fixation forceps  
1x2 straight, 0.50 mm teeth  
**RB50-700008**



Scleral forceps, straight,  
1x2 straight 0.70 mm teeth, with tying platform  
**RB50-700009**


### Conjunctival and Fixation Forceps

Conjunctival and fixation, forceps  
1x2 straight, 0.50 mm teeth  
**RB50-700010**



### Suturing forceps


Suturing forceps, straight,  
1x2 oblique teeth, with tying platform  
**RB50-700011-1** 0.12mm  
**RB50-700011-2** 0.35mm  
**RB50-700011-3** 0.60mm  
**RB50-700011-4** 0.90mm



## Forceps


### Jayle Corneal Tying Forceps

**RB50-5505** 0.15mm teeth  
**RB50-5506** 0.3mm teeth





### Jayle Corneal Tying Forceps

**RBT50-5507** 0.12mm teeth titanium  
**RBT50-5508** 0.3mm teeth titanium







### St.Martin Corneal Tying Forceps

**RB50-5509**  
 0.15mm teeth  
**RB50-5510**  
 0.3mm teeth  
**RBT50-5511** 0.15mm teeth titanium  
**RBT50-5512** 0.3mm teeth titanium



### Bonn Micro Tying Forceps

**RB50-5513**  
 corneal 0.12mm teeth  
**RB50-5515**  
 iris 0.12mm teeth  
**RBT50-5516** iris  
 0.12mm teeth titanium



### Mcpherson

Tissue Forceps


1x2 teeth, 0.4mm with tying platform

**RB50-2300** straight

**RB50-2310** angled

8mm long 1x2 teeth, 0.3mm

**RB50-2330** angled




### Hoffer - Mcpherson

Forceps very delicate

11mm long jaws

1x2 teeth, 0.12mm

**RB50-2350**



### Castroviejo Suturing Forceps

1x2 teeth with tying platform

**RB50-2500** 0.12mm

**RB50-2510** 0.3mm

**RB50-2520** 0.5mm

**RB50-2530** 0.9mm

**RBT50-2500** 0.12mm, titanium



### Castroviejo Suturing Forceps

1x2 teeth with tying platform


**RB50-25500** 0.12mm

**RB50-25501** 0.3mm

**RB50-25502** 0.5mm

**RB50-25503** 0.9mm

**RBT50-25504** 0.12mm, titanium



### Castroviejo Suturing Forceps


1x2 teeth with tying platform

**RB50-25300** 0.12mm

**RB50-25300-1** 0.3mm

**RB50-25300-2** 0.5mm


**RB50-25300-3** 0.9mm



## Forceps

### Moody

Fixation Forceps  
1x2 teeth, 0.5mm  
with lock  
**RB50-2553** left  
**RB50-2554** right




### Castroviejo-Colibri

Castroviejo-Colibri  
type corneal forceps 0.12mm, 1x2  
**RB50-1930**




### Girard

Corneoscleral Forceps  
1x2 teeth, 0.12mm  
with tying platform  
**RB50-2900** right hand  
**RB50-2910** left hand




### Stern-Castroviejo

Fixation Forceps  
1x2 teeth, 0.5mm with lock  
**RB50-2552**




### Bores

Corneal Fixation Forceps  
U - Shaped, 3mm spread  
1x2 teeth at each tip, 0.12mm  
**RB50-3250**




### Bishop-Harmon

Tissue Forceps  
1x2 teeth  
**RB50-3500** delicate  
**RB50-3510** standard  
**RB50-3520** heavy




### Bishop-Harmon

Dressing Forceps  
**RB50-3550** criss-cross serrated tips  
**RB50-3550-1** serrated 0.5  
**RB50-3550-2** serrated 0.8



### Fixation Forceps

1x2 teeth, 0.3mm  
**RB50-3560** sauer  
**RB50-3562** manhattan Eye & Ear  
with tying platform




### Dressing Forceps

Serrated, delicate  
**RB50-4000** straight  
**RB50-4010** curved



### Tissue Forceps


1x2 teeth  
**RB50-4100** delicate, straight  
**RB50-4110** delicate, curved  
**RB50-4140** standard, straight  
**RB50-4141** standard, curved



## Forceps

### Adson

straight, 5" (125mm) long  
**RB50-4150** serrated jaws, delicate  
**RB50-4152** serrated jaws, standard  
**RB50-4170** 1x2 teeth, delicate  
**RB50-4172** 1x2 teeth, standard  
**RBT50-4152** serrated titanium  
**RBT50-4170** 1x2 delicate titanium



### Bayonet

Dressing Forceps serrated jaws 5  
 3/4" (146mm) long  
**RB50-4180**




### Lester

Fixation Forceps straight  
**RB50-4300** 1x2 teeth  
**RB50-4320** 2x3 teeth



### Fixation Forceps

1x2 teeth  
**RB50-4500** Bracken delicate teeth  
**RB50-4550** Elschnig, standard teeth



### Troutman

Superior Rectus  
 Forceps 1x2 teeth angled  
**RB50-4700**



### Tying Forceps

For 8-0 to 11-0 sutures forceps

#### McPherson

**RB50-5000** 5mm long tying straight

**RBT50-5000** 5mm long tying straight, titanium

**RB50-5010** 5mm long tying angled

**RBT50-5010** angled, 5mm long tying titanium

#### Kelman - McPherson forceps

**RB50-5030** 7.5mm long tying angled

**RBT50-5030** 7.5mm long tying angled, titanium


#### Bechert Forceps

**RB50-5049** 10mm long tying angled


#### Gaskin Forceps

**RB50-5050** 11mm long tying angled

**RBT50-5050** 11mm long tying angled Titanium



RB50-5000 RB50-5010



RB50-5030 RB50-5049 RB50-5050



### Haberle-McPerson


Lens Holding Forceps

7.5mm long jaws

with groove, angled

For grasping IOL haptic

**RB50-5040**



### Kurstin

Flap Stretching Forceps

disc shaped

serrated jaws designed for

stretching the flap to

eliminate striae.


**RB50-5056**

LASIK Flap Forceps

disc shaped

serrated jaws

**RB50-5058**



RB50-5056 RB50-5058






## Forceps

### Inamura

Capsulorhexis Forceps delicate tapered shanks to seal the incision during continuous curvilinear capsulorhexis Excellent for MICS  
**RB50-5070**



### Utrata

Capsulorhexis Forceps very delicate, short handle triangular grasping tips extremely thin with landmarks at 2.5mm and 5.0mm  
**RB50-5071** straight




### Masket

Capsulorhexis Forceps very delicate, short handle triangular grasping tips extremely thin with landmarks at 2.5mm and 5.0mm  
**RB50-5074** vaulted



### Kraff

Capsule Tag Forceps  
 11mm long shanks  
 1.5mm long grasping tips  
**RB50-5080**



### Utrata Capsulorhexis Forceps

Sharp Tips, 80mm  
**RB50-2002S**



**Utrata**

Capsulorhexis Forceps  
 very delicate  
 triangular grasping tips  
 extremely thin  
 11mm long shanks  
**RB50-5081** stainless steel  
**RBT50-5081** titanium




**Utrata**

Capsulorhexis Forceps very  
 delicate  
 triangular grasping tips  
 extremely thin  
 11mm long shanks  
 round handle with guide pin  
 for alignment  
**RB50-5082**



**Nevyas**

Capsulorhexis Forceps  
 cystotome shaped tips  
 extremely thin  
 11mm long shanks  
 For tearing and grasping  
 the anterior capsule.  
**RB50-5083** stainless steel  
**RBT50-5083** titanium




**Masket**

Capsulorhexis Forceps  
 very delicate  
 triangular grasping tips  
 extremely thin 11mm long  
 vaulted shanks  
**RB50-5084**  
**RBT50-5084** titanium



**Inamura**

Capsulorhexis Forceps  
 round handle cross action  
 2 mm incision  
**RB50-50070**



## Forceps

### Masket

Capsulorhexis Forceps  
very delicate triangular  
grasping tips extremely  
thin 11mm long vaulted  
shanks round handle with  
guide pin for alignment

**RB50-5086**


**RBT50-5086** titanium



### Capsulorhexis Forceps

Corydon forceps  
ergonomic curved, cystotome tips,  
11.5 cm long

**RB50-17163**



### Capsulorhexis Forceps

Capsulorhexis Forceps  
angled, 12 mm graduated slightly curved  
arms, sharp tips, ergonomic handle

**RB50-600002** 11mm angled


**RB50-600003** 12mm angled



### Capsulorhexis Forceps

Round handle  
11mm long shanks straight


**RB50-405081**



### Capsulorhexis Forceps

Round handle, curved,  
cystotome tips, 11.5cm long

**RB50-500091**




### Inamura Capsulorhexis Forceps

Sharp pointed serrated tips angled 90° from shaft  
 Curved shaft, tip to pivot point 8.5mm  
 Cross action tips, 1.2mm width at pivot box  
 Round handle, overall length 111mm desired size of capsulorhexis

**RBT50-00786**

**RBT50-00786-1** marks on shaft at 2.5mm and 5.0mm denote



### Capsulorhexis forceps

Capsulorhexis forceps, lightweight, angled, 12 mm, slightly curved arms, sharp tips,

**RBT50-80001234** titanium

Capsulorhexis forceps, lightweight, angled, 12 mm, slightly curved arms, strongly angled, sharp tips,


**RBT50-80001234-1** titanium

Capsulorhexis forceps, lightweight, angled, 11 mm, slightly curved arms, notched tips

**RBT50-80001234-2** titanium

Capsulorhexis forceps, lightweight, angled, 11 mm, slightly curved arms, sharp, strongly angled tips

**RBT50-80001234-3** titanium



### Capsulorhexis forceps

Capsulorhexis forceps, angled, 11 mm slightly curved arms, sharp, notched tips

**RB50-700003**



### Capsulorhexis forceps

Capsulorhexis forceps, angled, 11 mm slightly curved arms, sharp tips ergonomic handle

**RBT50-700002** titanium



## Forceps

### Capsulorhexis forceps

Capsulorhexis forceps  
flat handle, 11 mm sharp tips,  
micro incision

**RBT50-700001** titanium



### MICS capsulorhexis forceps

MICS capsulorhexis forceps,  
controlled opening for 2 mm incisions

**RBT50-700018** titanium



### Capsulorhexis forceps for micro coaxial phaco

Capsulorhexis forceps with 2mm  
controlled opening for 2mm incisions

**RBT50-500099** titanium



### MICS/MCCS Capsulorhexis Forceps

MICS/MCCS Capsulorhexis  
with 2 mm controlled opening,  
for 2 mm incisions


**RB50-600000**



### Capsulorhexis Forceps

extremely delicate

**RB50-31300**




### Prechopper

Combo Prechopper,  
for grades 1-4 nuclei,  
2 mm opening, sharp top side  
**RB50-60005-1**




### Prechopper


Akahoshi Super Combo III Prechopper  
For nucleus grades 1 to 3.  
is suitable for the vertical Karate  
Prechop technique within a 2mm incision.  
Prechopping can be completed without  
nuclear support for nuclei grade 1 to 3  
**RB50-7233**



Akahoshi Hybrid Combo prechopper  
For Grades 1-4 Nuclei  
**RB50-7235**




Akahoshi Universal Nucleus Splitter  
**RB50-7236**



### Akahoshi

Phaco PreChopper  
4.5mm long delicate jaws  
with sharp tips  
For cracking, dividing and rotating  
the nucleus without sculpting  
**RB50-7230**



### Akahoshi


Combo PreChopper  
For cracking, dividing and rotating  
both hard and soft nuclei without  
sculpting  
**RB50-7232**



## Forceps


### Anis

Tying Forceps  
extra delicate  
smooth jaws  
For 9-0 to 11-0 sutures.  
**RB50-5112** straight  
**RB50-5113** curved




### Freeman - Jaffe

Lens Holding Forceps  
with longitudinal groove  
For grasping IOL haptic.  
**RB50-5190**





### Jaffe

Tying Forceps  
very delicate  
smooth jaws  
6mm long platform  
For 8-0 to 11-0 sutures.  
**RB50-5200** straight  
**RB50-5210** curved  
**RBT50-5200** straight titanium  
**RBT50-5210** curved titanium





### Gerl

Tying Forceps  
curved shanks  
6mm long platform  
For 8-0 to 11-0 sutures.  
**RBT50-5215** titanium



### Tennant

Tying Forceps  
extra delicate smooth jaws  
with 6mm long platform  
round handle with guide pin  
For 9-0 to 11-0 sutures.  
**RB50-5220** straight  
**RBT50-5220** straight, titanium  
**RB50-5230** curved  
**RBT50-5230** curved, titanium



### Maumenee

Tying Forceps  
concave-convex jaws  
round handle with guide pin  
For 8-0 to 11-0 sutures.

**RBT50-5234** straight titanium

**RBT50-5235** angled titanium

**RB50-5234** straight

**RB50-5235** angled



Jaws Cross Section



### Harms


Tying Forceps  
standard, smooth jaws  
for 7-0 to 9-0 sutures.

**RB50-5400** straight

**RB50-5410** curved

**RBT50-5400** straight titanium

**RBT50-5410** curved titanium




### Tuebingen Model

Tying forceps straight

**RB50-19030** S.S

**RBT50-19030** titanium




### Tuebingen Model

Tying forceps curved

**RB50-19040** S.S

**RBT50-19040** titanium



### McPherson Tying Forceps


Recommended for 8-0 to 10-0 sutures

**RB50-700019** 4mm tying platform straight

**RB50-700019-1** 4mm tying platform curved

**RB50-700019-2** 4mm tying platform angled

**RB50-700019-3** 8mm tying platform angled






## Forceps

### Kelmann forceps

Kelmann forceps angled ergonomic handle


**RB50-600007** 11 mm tying platform  
**RB50-600007-1** 7 mm tying platform



### Troutmann forceps

Troutmann forceps very delicate ergonomic handle

**RB50-600075** 7 mm tying platform str  
**RB50-600008** 7 mm tying platform cvd



### Troutman forceps

Straight, delicate, bullet shape.  
 Curved, delicate, bullet shape.


**RB50-9604**



### IOL Loading forceps

Loading forceps for foldable lenses into the cartridge, 11 mm angle to tip, cylindrical arms


**RBT50-700014** titanium



### Golf Ball Tying Forceps


Golf ball tying forceps titanium

**RBT50-700016-1** angled 8mm  
**RBT50-700016-2** straight 5mm  
**RBT50-700016-3** curved 5mm  
**RBT50-700016-4** delicate str 5mm





### Tying Forceps

Troutman forceps 6mm  
**RB50-55410** straight  
**RB50-55420** curved




### Kelman McPherson

French pattern tying forceps  
11mm angled  
**RB50-500088**




### McPherson

**RB50-511188** Tying forceps 5mm str  
**RB50-500089** Tying forceps 5mm angled  
**RB50-500090** Tying forceps 10mm angled




### Fixation forceps

Fixation forceps, straight,  
**RB50-700012-1** 3mm  
**RB50-700012-2** 4mm



### Green Fixation Forceps


Green fixation forceps  
**RB50-700013**



## Forceps


### Helveston

Tying Forceps  
standard, smooth jaws  
For 4-0 to 6-0 sutures.  
**RB50-5420**



### Fine

Tying Forceps  
standard, smooth jaws  
For 7-0 to 9-0 sutures.  
**RB50-5550**




### Arita

Meibomian Gland Compressor  
broad, smoothly polished jaws  
**RB50-5900**



### Cilia Forceps with spring tension

**RB50-6205**  
**RBT50-6205 titanium**



### Barraquer

Cilia Forceps lightly angled  
Textured jaws provide  
superior grasping capabilities.  
**RB50-6000**



### Barraquer Cilia Forceps

**RB50-6001**



**Beer Cilia Forceps**  
RB50-9300



**Cilia Forceps**

Cilia Forceps  
slanted jaws gold plated  
RB50-6020




**Gradle**

Cilia Forceps smooth jaws  
RB50-6100




**Bonaccolto**

Conjunctiva Forceps  
longitudinally serrated jaws  
with cross serrations at tips  
RB50-6200 1.2mm wide  
RB50-6210 1.7mm wide



**Giannetti**


MICS Capsulorhexis Forceps  
ultra delicate shanks and tips  
maximum jaw opening 1.75mm  
For performing capsulorhexis  
through a micro incision.  
RB50-5090



## Forceps


### Moorfields Conjunctival Forceps

- RB50-5500F standard
- RB50-5501F fine
- RBT50-5500F standard titanium
- RBT50-5501F fine titanium




### Stallard Tying Forceps

RB50-5502



### Rycroft Corneal Tying Forceps

- RB50-5503 plain
- RB50-5504 0.3mm teeth



### Scleral Plug Forceps

RB50-5517




### Vessel Dilator

11 cm long, flat handle 9 mm wide,  
angulated 10°  
**RB50-6590st**



### Micro Forcep

15 cm long, round handle  
**RB50-6591st** 6mm tying platform 1x2, 0.30 straight  
**RB50-6592st** 6mm tying platform straight  
**RB50-6593st** 6mm tying platform curved  
**RB50-6594st** vessel dilator



### Nugent

Utility Forceps  
**RB50-6214** smooth jaws  
**RB50-6215** serrated jaws



### Foreign

Foreign Body Forceps straight  
**RB50-6400**



### Utility & Splinter

Forceps pointed, short jaws  
Suitable for epilation.  
**RB50-6500 #H**



## Forceps

### Jewelers

Forceps  
standard, straight  
**RB50-6510 #1**



### Jewelers

Forceps  
delicate, straight  
**RB50-6520 #3-C**



### Jewelers

Forceps  
delicate, straight  
**RB50-6530 #3**



### Jewelers

Forceps  
delicate, straight  
**RB50-6540 #4**



### Jewelers

Forceps  
extra delicate, straight  
**RB50-6550 #5**




**Osher**

Iris Forceps  
very fine textured jaws  
angled  
**RB50-6558**



**Jewelers**

Forceps  
extra delicate, angled  
**RB50-6560 #5/45**




**Jewelers**

Forceps  
curved  
**RB50-6570 #7**



**Suture**

Removal Forceps  
very fine pointed tips with textured  
jaws Ideal for grasping wet sutures  
**RB50-6575**



**Jaffe**

Utility Forceps  
very delicate smooth jaws curved  
**RB50-6600**





## Forceps

### Arruga

Capsule Forceps  
RB50-7000



### Rigid IOL implantation forceps

Rigid IOL implantation forceps, angled  
RB50-700015



### Implantation forceps

Implantation forceps  
for 5.5 and 6.0 mm acrylic IOLs  
RB50-600006



### Ernest

Nucleus Cracker  
extra delicate cross-action  
serrated paddle tips  
For cracking the nucleus  
through a 2.5mm incision  
RB50-7240




### Ernest

Micro Nucleus Cracker  
extra delicate cross-action  
serrated paddle tips  
For cracking the nucleus  
through a 2mm incision.  
RB50-7242



**Dodick**

Nucleus Cracker  
cross-action  
longitudinally serrated  
paddle shaped jaws  
For cracking the nucleus  
through a 3mm incision  
**RB50-7244**



**Hu Femto Forceps**

Blunted 0.3mm tips  
cross-action mechanism  
For opening incisions  
created with a  
femtosecond laser  
**RB50-7245**



**Alfonso**

Nucleus Grasping Forceps  
two short rows of  
delicate interlocking teeth  
For grasping and extracting  
nucleus fragments  
**RB50-7253**



**Coaxial forceps**

Coaxial forceps for endothelial  
lenticule insertion in DSEK technique,  
atraumatic  
**RBT50-3839.5**



**Corneal-scleral punch**


Corneal-scleral punch,  
straight, Ø 1.5 mm Ø 1.7 mm  
**RB50-2865**



## Forceps


### Stripping forceps

Stripping forceps for the Schlemm's canal, straight  
**RB50-3276.7**




### Stripping forceps

Stripping forceps for the Schlemm canal, curved tip  
**RB50-3279**






### Dodick

IOL folding forceps  
 for soft IOLs  
**RB50-8210**



### Livernois

IOL folding forceps  
 concave jaw surfaces for soft IOLs  
**RB50-8215**



**Faulkner**

Lens Holding Forceps  
 round, smooth jaws  
 For holding soft IOLs  
**RB50-8230**




**Buratto**

Soft IOL Inserting Forceps  
 extra delicate highly polished  
 convex jaws for foldable acrylic and  
 silicone IOLs  
**RB50-8268**



**Spaleck**

Soft IOL Inserting Forceps  
 delicate highly polished convex jaws  
 For foldable acrylic and silicone IOLs  
**RB50-8270**



**Ernest-McDonald**

IOL Inserting Forceps  
 polished bi-convex jaws  
 patented cross-action for  
 acrylic IOLs  
**RB50-8272**



JAWS CROSS SECTION



**Bechert**

Lens Holding Forceps  
 very delicate flat jaws  
 gently curved  
**RB50-8476**



## Forceps

### Blaydes

Lens Holding Forceps  
very delicate narrow jaws  
angled 45°  
**RB50-8500**




### Shepard

Lens Holding Forceps  
duckbill shaped jaws  
gently curved  
**RB50-8600**



### Clayman

Lens Holding Forceps  
delicate angled jaws  
gently curved  
**RB50-8700**



### Jameson

Muscle Forceps  
with slide lock  
**RB50-9000** 4 teeth, child size, right  
**RB50-9010** 4 teeth, child size, left  
**RB50-9020** 6 teeth, adult size, right  
**RB50-9030** 6 teeth, adult size, left



### Berke

Ptosis Forceps  
with slide lock  
longitudinally grooved jaws  
**RB50-9300** 20mm long  
**RB50-9310** 27mm long




### Steinert

IOL Inserting Forceps  
conically shaped jaws  
with flat inside surfaces  
Used for folding and  
inserting soft IOLs.  
**RB50-8219**




### Livernois

IOL Holding Forceps  
highly polished round jaws  
For holding a soft IOL  
while folding with an  
inserting forceps.  
**RB50-8220** straight  
**RB50-8221** angled




### Livernois-McDonald

IOL Inserting Forceps  
5mm long concave jaws  
patented cross-action  
For inserting silicone IOLs  
**RB50-8225**




### Ernest-McDonald

IOL Inserting Forceps  
5mm long straight jaws  
low profile  
patented cross-action  
For inserting silicone IOLs  
**RB50-8228**



### Ernest-McDonald II

Lens Cartridge Loading Forceps  
straight, 20 gauge jaws feature a  
series of concave and convex surfaces  
**RB50-8229**




## Forceps

### Chalazion Forceps

round

**RB50-9450** Lambert, 8mm I.D.

**RB50-9500** Hunt, 12mm I.D.




### Chalazion Forceps

oval

**RB50-9550** Baird, 11 x 8mm I.D. small

**RB50-9600** Francis, 14 x 11mm I.D. medium

**RB50-9650** Desmarres, 20 x 13mm I.D. large



### Desmarres Forceps

Oval

**RB50-2125** 25mm



### Watzke

Sleeve Spreading Forceps

**RB50-7950**




**Hartmann**

Mosquito Forceps  
extra delicate  
3 3/4" (95mm) long  
**RB50-9900** straight  
**RB50-9910** curved




**Halsted**

Mosquito Forceps  
delicate  
5" (125mm) long  
**RB50-9940** straight  
**RB50-9950** curved



**Halsted**

Mosquito Forceps  
delicate 1x2  
5" (125mm) long  
**RB50-9970** straight  
**RB50-9980** curved





## Forceps

### Backhaus

Towel Clamp

**RB50-9964** 3½" (90mm)

**RB50-9965** 5¼" (135mm)

**RBT50-9964** 3½" (90mm) titanium

**RBT50-9965** 5¼" (135mm) titanium



### Baby Jones

Towel Clamp

Cross Action

**RB50-9960** (55mm) long



### Schaedel

Towel Clamp

Cross Action

**RB50-9961** (90mm) long



### Serrefine

serrated jaws

**RB50-9850** small (38mm)



### Dieffenbach Serrefine

**RB50-9851** Straight, 50 mm

**RB50-9852** curved, 50 mm




### Pean Forceps

Mosquito forceps,  
straight, serrated jaws  
**RB50-2720**



### Ramsey Forceps

Ramsey sponge forceps  
**RB50-9920** 18cm  
**RB50-9921** 25cm



### Frazier

Suction and wash tube,  
with hole on/off, 12cm long  
**RB20-1861**

Suction and wash tube,  
with hole on/off, 17cm long  
**RB20-1862**



rbmedical

# RB60

• Needle Holders

## Needle Holder

### McPherson

Needle Holder  
very delicate jaws  
**RB60-1000** with lock straight  
**RB60-1010** without lock straight  
**RB60-1020** with lock curved  
**RB60-1030** without lock curved



### Jaffe

Needle Holder  
delicate conical jaws  
**RB60-1100** with lock straight  
**RB60-1120** with lock curved  
**RB60-1130** without lock curved




### Castroviejo

Needle Holder  
delicate jaws straight  
**RB60-1500** with lock  
**RB60-1510** without lock



### Castroviejo

**RB60-1600** with lock  
**RB60-1610** without lock



### Silcock Needle Holder,catch

**RB60-1611**



## Needle Holder

### Castroviejo

Needle Holder  
delicate jaws  
curved  
**RB60-1520** with lock  
**RB60-1530** without lock



### Castroviejo

Needle Holder  
standard jaws  
straight, with lock  
**RB60-1560**  
**RBTC60-1560** Tungsten Tip TC



### Castroviejo

Needle Holder  
standard jaws curved  
**RB60-1580** with lock  
**RB60-1581** without lock  
**RBT60-1580** with lock  
**RBT60-1581** without lock




### Kalt

Needle Holder  
delicate smooth jaws  
straight, with lock  
**RB60-1800**



### Drews

Needle Holder  
extra delicate jaws  
**RB60-2990** straight, with lock  
**RB60-2991** straight, without lock  
**RB60-2992** curved, with lock  
**RB60-2993** curved, without lock



### Castroviejo

Needle Holder  
delicate jaws titanium  
**RBT60-1500** with lock str  
**RBT60-1520** with lock cvd  
**RBT60-1530** without lock cvd  
**RBT60-1560** with lock standard



### Barraquer

Needle Holder  
short model  
very delicate jaws  
curved  
**RB60-3020** with lock  
**RB60-3030** without lock



### Troutman

Needle Holder  
delicate jaws  
strongly curved  
**RB60-3250** without lock  
**RB60-3260** with lock



### Barraquer

Needle Holder  
very delicate jaws  
curved  
**RB60-3300** with lock  
**RB60-3310** without lock  
**RBT60-3300** titanium, with lock  
**RBT60-3310** titanium, without lock



### Barraquer

Needle Holder  
delicate jaws  
straight  
**RB60-3500** with lock  
**RB60-3510** without lock



## Needle Holder

### Barraquer

Needle Holder  
delicate jaws  
curved  
**RB60-3520** with lock  
**RB60-3530** without lock

Needle Holder  
delicate jaws  
curved large  
**RB60-3540** with lock  
**RB60-3550** without lock



### Barraquer

Needle Holder  
standard jaws  
curved  
**RB60-3620** with lock  
**RB60-3630** without lock  
**RBTC60-3620** Tungsten Carbide with lock  
**RBTC60-3630** Tungsten Carbide without lock  
**RBT60-3630** with lock titanium  
**RBT60-3630** without lock titanium




### Anis

Needle Holder  
extra delicate jaws  
**RB60-3750** with lock straight  
**RB60-3751** without lock straight  
**RB60-3752** with lock curved  
**RB60-3753** without lock curved



### Micro Needle Holder

Needle Holder  
extra delicate jaws  
straight  
**RB60-3800** with lock  
**RB60-3810** without lock  
**RBT60-3800** titanium, with lock  
**RBT60-3810** titanium, without lock




### Barraquer ( Golf Pattern )

Needle Holder  
short model  
very delicate jaws  
Curved  
**RB60-1620** with lock  
**RB60-1630** without lock



### Barraquer

Needle Holder  
short model  
very delicate jaws  
curved  
**RB60-1640** with lock  
**RB60-1650** without lock



### Micro Needle Holder

Needle Holder  
delicate jaws straight  
**RB60-1660** without lock  
**RB60-1670** with lock



### Needle Holder

delicate,  
curved, without lock,  
overall length: 115 mm, titanium  
Ideal for 9-0 to 10-0 sutures  
**RB60-1680** with lock  
**RB60-1690** without lock  
**RBT60-1680** titanium, with lock  
**RBT60-1690** titanium, without lock



### Barraquer

Needle Holder  
delicate jaws  
straight  
**RB60-1700** with lock  
**RB60-1710** without lock  
**RBT60-1700** with lock, titanium  
**RBT60-1710** without lock, titanium





# RB70

- Lacrimal Instruments
- Bimanual I/A
- Quickswitch I/A System
- Infusion Handles
- Rongeurs
- Nasal Speculums
- Elevators
- Chisels
- Mallets
- Probes

**Kerrison**

Bone Nibbling Rongeur  
9mm opening

**RB70-1800** size 0, 3.0mm wide

**RB70-1801** size 1, 3.5mm wide

**RB70-1802** size 2, 4.0mm wide

**RB70-1803** size 3, 4.5mm wide



2/3

**Citelli**

Bone Nibbling Rongeur  
7mm opening

2.0mm wide

**RB70-1805**



2/3

**Beyer**

Rongeur  
single action

**RB70-1810**



2/3

## Lacrimal Instruments

### Nasal

Speculum  
RB70-1820 infant  
RB70-1821 child  
RB70-1822 adult



2/3

### Freer

Periosteal Elevator  
one end blunt  
other end semi-sharp  
4.5mm wide blades  
RB70-1830



2/3

### Joseph

Periosteal Elevator  
gently curved tip  
RB70-1835



3/4

### West

Lacrimal Chisel  
straight  
RB70-1840



3/4

### Freer

Lacrimal Chisel  
curved  
RB70-1845



**Mallet**


lead filled  
stainless steel head  
**RB70-1860**



2/3

**Harms**

Trabeculotomy Probe  
**RB70-1000** right  
**RB70-1010** left



**Infant**

Lacrimal Dilator  
very fine blunt tip  
**RB70-2000**




**Muldoon**

Lacrimal Dilator  
**RB70-2100**



**Wilder**

Lacrimal Dilator  
**RB70-2110** size 1, fine taper  
**RB70-2120** size 2, medium taper  
**RB70-2130** size 3, heavy taper



## Lacrimal Instruments


### Castroviejo

Lacrimal Dilator  
double-end  
fine and  
medium taper  
**RB70-2300**



### Bowman

Lacrimal Probe  
malleable stainless steel  
**RB70-2500** set of 6  
**RB70-2510** size 0000 & 000  
**RB70-2520** size 00 & 0  
**RB70-2530** size 1 & 2  
**RB70-2540** size 3 & 4  
**RB70-2550** size 5 & 6  
**RB70-2560** size 7 & 8



3/4

### Catalano

Lacrimal Intubation Set  
includes 2 blunt  
0.7mm diameter probes  
attached to 30.5cm long  
silicone tubing  
.012" I.D., .025" O.D.  
**RB70-2610** set  
**RB70-2602** silicone tubing only



3/4

### Crawford

Lacrimal Intubation Set  
includes 2 olive tip  
0.6mm diameter probes  
attached to 30.5cm long  
silicone tubing  
.012" I.D., .025" O.D.  
**RB70-2620** set  
**RB70-2602** silicone tubing only



3/4

**Crawford**

Retrieving Hook  
for Crawford lacrimal  
intubation probes  
**RB70-2623**



**Pigtail**

Probe  
with suture holes  
**RB70-2900**



**Rollet Rugine**

**RB70-2910**



**Nettleship Lacrimal Dilator**

Nettleship lacrimal dilator  
single end  
**RB70-2920**  
Nettleship lacrimal dilator  
double end  
**RB70-2930**



**Howarth Periosteal Elevator**

Howarth Periosteal Elevator,  
for DCR, 215mm length  
**RB70-2940**




**Fison Retractor**

Fison Retractor, for detachment  
**RB70-2950**




# Cannulas



## Lacrimal

Cannula  
malleable tip  
23 gauge  
**RB7-3000** straight  
**RB7-3010** curved



## Air Injection

Cannula  
angled at 7mm  
**RB70-3502** 30 gauge  
**RB70-3522** 27 gauge



Bottom View



## Dalk

Air Injection Cannula  
flattened, 27 gauge  
vaulted 0.2mm  
posterior port near tip  
**RB70-3590** 27 gauge



Top View



Bottom View



## Dalk

DALK Cannula flattened  
"cobra-shaped" tip  
0.2mm posterior port  
**RB70-3595** 27 gauge



## Lasek

Alcohol Cannula  
delicate olive shaped tip  
angled at 10mm  
25 gauge  
**RB70-3805**





## Stevens

Sub-Tenon's Anesthesia  
Cannula gently curved,  
flat tip 19 gauge  
**RB70-4005**



## Gess

Irrigating-Aspirating  
Cannula U-shaped,  
with end opening  
aspiration port,  
23 gauge Supplied  
with silicone tubing  
and adaptor  
**RB70-4304** curved right  
**RB70-4305** curved left



## Simcoe


Irrigating-Aspirating  
Cannula original model,  
15mm long,  
0.3mm diameter  
aspiration port  
**RB70-4300**



## Corbin


Sub-Tenon's Anesthesia  
Cannula smooth rounded  
tip 0.3mm side port  
23 gauge  
**RB70-4008**

Cannulas




**Reverse Simcoe**

Irrigating-Aspirating  
Cannula 15mm long  
0.3mm diameter  
aspiration port front  
opening irrigation port  
23 gauge  
**RB70-4302**




**Simcoe**

Irrigating-Aspirating Cannula  
U-shaped, with 0.3mm anterior  
aspiration port, 23 gauge  
**RB70-4308** curved right  
**RB70-4309** curved left




**Strauss-Simcoe**

MSICS I/A Cannula  
23 gauge, 0.3mm port  
elongated and gently curved  
with 15° angled tip.  
**RB70-4320**




**Kelman**

Cystotome  
sharp point and  
sharp cutting edge  
22 gauge  
**RB70-4600**




**McIntyre**

Cystotome  
angled at 10mm  
sharp point  
23 gauge  
**RB70-4680**




**Lieppman**

Cystotome  
extra small sharp point  
and sharp cutting edge  
22 gauge  
**RB70-4700**




**Shammas**

Posterior Capsule  
Polishing Curette  
1mm diameter cup  
with textured base  
22 gauge  
**RB70-4850**




**Shapiro**

Posterior Capsule  
Polisher carbide  
impregnated tip dual  
0.3mm side ports  
curved, 23 gauge  
**RB70-4906**



**Kratz**

Cystotome  
sharp point and  
sharp cutting edge  
formed, 22 gauge  
**RB70-4601**



**Jensen**

Posterior Capsule  
Polisher carbide  
impregnated on posterior  
surface of olive tip  
front opening  
**RB70-4920** 25 gauge  
**RB70-4930** 23 gauge



## Cannulas

### Bishop-Harmon

Anterior Chamber Irrigator includes 19 gauge angled cannula, adaptor and blue silicone bulb

- RB70-5000** complete
- RB70-5005** silicone bulb and adaptor
- RB70-5010** silicone bulb
- RB70-5020** adaptor
- RB70-5050** 19 gauge cannula
- RB70-5052** 21 gauge cannula



### Lu

Posterior Capsule Polisher flattened olive shaped tip carbide impregnated on the posterior surface front opening, 23 gauge  
**RB70-4933**



### Dahan

Nucleus Prolapsing Cannula angled at 4mm miniature olive tip with 0.15mm end opening 25 gauge  
**RB70-5058**



### Hydrodissection

Cannula angled at 8mm flat tip, front opening 25 gauge  
**RB70-5054**




### Irrigating Cannula

dual 0.4mm side ports 23 gauge thin wall  
**RB70-5095**



### Fukasaku

LASIK Cannula low profile tip, end opening 23 gauge  
**RB70-5108**



### Banaji

LASIK Cannula bullet shaped tip 6 ports for multi-directional irrigation 25 gauge  
**RB70-5112** smooth  
**RB70-5113** textured




### Buratto

LASIK Cannula slightly flattened one front and two side ports  
**RB70-5117** 25 gauge  
**RB70-5118** 23 gauge



### McIntyre

Anterior Chamber Cannula smooth blunt tip, front opening  
**RB70-5120** 23 gauge, straight  
**RB70-5130** 23 gauge, angled  
**RB70-5140** 26 gauge, straight  
**RB70-5150** 26 gauge, angled  
**RB70-5160** 30 gauge, straight  
**RB70-5170** 30 gauge, angled



### Gimbel

Cortex Aspiring Cannula end opening 30 gauge  
**RB70-5174**




**Vidaurri**

LASIK Cannula  
double-armed  
8 irrigating ports  
25 gauge  
**RB70-5180** vertical  
**RB70-5181** horizontal



**Simcoe**

Cortex Extractor  
0.3mm side port  
13mm long, 23 gauge  
**RB70-5220** smooth tip  
**RB70-5221** textured tip



**Simcoe**

Cortex Extractor  
0.3mm side port  
21mm long, 23 gauge  
**RB70-5230** smooth tip  
**RB70-5231** textured tip




**McIntyre-Binkhorst**

Cannula, "J" shape  
smooth blunt tip  
end opening, 26 gauge  
**RB70-5440** angled right  
**RB70-5445** straight  
**RB70-5450** angled left




**Chang**

Hydrodissection Cannula  
flat tip angled 90°  
with beveled opening  
27 gauge  
**RB70-5464** standard tip  
**RB70-5466** long tip




**Knolle-Pearce**

Irrigating Vectis  
fluid flow through  
three front opening  
ports 23 gauge loop  
**RB70-5500**




**Simcoe**

Irrigating Lens Nucleus Loop  
two front opening ports  
finely serrated, strong curve  
25 gauge  
**RB70-5530** right  
**RB70-5531** left



**Simcoe**

Cortex Extractor  
0.4mm side port  
21mm long  
23 gauge thin wall  
**RB70-5232**



**Alfonso-McIntyre**

Irrigating Nucleus Spoon  
single 0.3mm  
front opening port  
**RB70-5532** 2.5mm wide  
**RB70-5533** 4.0mm wide

## I/A Bimanual Handles

### Infusion Handle


Infusion Handle  
 RB70-6805 male/male connectors  
 RB70-6806 male/female connectors



### Aspiration Handpieces

with male connector  
 side port facing into curve  
 magenta color titanium  
 titanium

- RBT70-5804** 23 gauge, 0.3mm port
- RBT70-5805** 23 gauge, 0.3mm port, textured tip
- RBT70-5810** 21 gauge, 0.35mm port
- RBT70-5811** 21 gauge, 0.35mm port, textured tip
- RBT70-5814** 21 gauge, 0.4mm port
- RBT70-5815** 21 gauge, 0.4mm port, textured tip
- RBT70-5816** 22 gauge, 0.3mm port, ova
- RBT70-5817** 22 gauge, 0.3mm port, oval, textured tip
- RBT70-5822** 20 gauge, 0.3mm port, oval



### Irrigation Handpieces

with female connector, blue color


- RBT70-5836** 23 gauge, dual 0.3mm side ports
- RBT70-5838** 22 gauge, dual 0.4mm side ports, oval
- RBT70-5840** 21 gauge, dual 0.5mm side ports

### Aspiration Handpieces Ø21G

Curved sandblasted tip  
 titanium handle  
**RBT70-2233-1**


### Irrigation Handpieces Ø21G

Curved tip, dual sideports,  
 titanium handle  
**RBT70-2232-2**



### Scoring Hook and Aspirating Handpiece

Scoring hook and aspirating handpiece  
 for DSEK procedure  
**RBT70-2183.14**



**Bipolar Forceps**

Jewelers Type  
pointed tips  
straight  
**RB80-7010** plain  
**RB80-7012** insulated



**Bipolar Forceps**

Bipolar Forceps  
**RB80-7020** Straight  
**RB80-7021** Angled



**Bipolar**

Coaptation Forceps  
pointed tips  
**RB80-7030**



**Bipolar**

Bayonet Forceps  
**RB80-7040**



**Bipolar Cord**

Autoclavable  
**RB80-7065**



## Sterilizing Cases & Trays

### Plastic Sterilizing Cases

Autoclavable plastic with silicone finger mat base, lid and mat are perforated for steam penetration


Catalog #	Size # of Instr	Sizes (mm)
<b>RB90-2020</b>	small 2-3	160 x 70 x 25



### Plastic Trays

Autoclavable plastic with silicone finger mat base, lid and mat are perforated for steam penetration

Catalog #	Size # of Instr	Sizes (mm)
<b>RB90-2030</b>	large 10-15	265 x 160 x 25
<b>RB90-2040</b>	large, double level 20-25	265 x 160 x 50



### Sterilization Perforated Case With Lid

S.S Sterilizing Case inox with silicone finger mat base, lid and mat are perforated for steam penetration

<b>RB90-22518356</b>	252 x 209 x 24mm silicone mat
<b>RB90-22518357</b>	189 x 209 x 24mm silicone mat
<b>RB90-22518358</b>	126 x 209 x 24mm silicone mat



## Sterilizing Cases & Tray

### Sterilization Tray with Lid

S.S Perforated stainless steel tray with lid with silicone fitting  
 Size: 240x250x38mm  
**RB90-22000**



### Fine Mesh Trays for Micro Instruments with Silicone Mat

These trays are especially designed for sterilizing delicate micro instruments.

The mesh size of these trays are 1.6mm, which helps to clean the instruments from all directions.

- RB90-4350** 160 X 70 X 20mm
- RB90-4351** 160 X 100 X 20mm
- RB90-4352** 200 X 125 X 20mm
- RB90-4353** 200 X 150 X 20mm
- RB90-4354** 265 X 160 X 20mm
- RB90-4355** 210 X 190 X 20mm
- RB90-4356** 250 X 250 X 20mm for Cataract Set



### Fine Mesh Cassettes With Silicone Fittings For Micro Instruments

- RB90-4248** 200 X 70 X 40 for 3 Micro Retinal Instruments
- RB90-4249** 200 X 90 X 40 for 5 Micro Instruments
- RB90-4250** 200 X 150 X 40 for 10 Micro Instruments
- RB90-4251** 200 X 210 X 40 for 15 Micro Instruments



## Sterilizing Cases & Tray

### Micro Instruments Tray

Tray with silicone finger mat base, lid and mat are perforated for steam penetration

<b>RB90-4243</b>	for 5 instr	185x60x23mm
<b>RB90-4244</b>	for 10 instr	185x30x23mm
<b>RB90-4245</b>	for 15 instr	185x205x45mm



### Ophthalmic tray

S.S Ophthalmic tray with silicone mat

<b>RB90-35000</b>	for 5 instruments
<b>RB90-35000-1</b>	for 10 instruments
<b>RB90-35000-2</b>	for 15 to 20 instruments



### Full Perforated Tray

Tray Specially designed for small parts and burrs

<b>RB90-4944</b>	40 X 40 X 30mm
<b>RB90-4945</b>	80 X 80 X 50mm
<b>RB90-4946</b>	90 X 90 X 50mm



# rbmedical

RB Medical Engineering Ltd  
Unit 2, Alton Road Industrial Estate  
Ross-on-Wye, Herefordshire  
HR9 5NS

T: +44 (0)1989 563958

F: +44 (0)1989 768267

email: [info@rbmedical.co.uk](mailto:info@rbmedical.co.uk)





