

# Speculum Insulated with smoke tube

Reusable



Product Code	Description - Insulated Cusco Speculum with Smoke Tube	Dimensions mm	
		A	B
DA160-V	Extra small	16	86
DA160-S	Small	25	95
DA160-SL	Small Long	25	115
DA160-SEL	Small Extra Long (Winterton size)	25	125
DA160-M	Medium	30	102
DA160-ML	Medium Long	30	115
DA160-MEL	Medium Extra Long	30	125
DA160-L	Large	35	115
DA160-LL	Large Long	35	125
DA160-EXLEXL	Extra Large Extra Long	38	135





Providing the instruments for good health

[www.rbmedical.co.uk](http://www.rbmedical.co.uk)

## **Key Benefits:**

- Extensive range of sizes available
- High quality construction and insulation for patient and user safety.
- Designed for improved patient comfort.
- Self retaining
- Convenient and easy to use
- Easy to clean and fully autoclavable

Our speculums and gynaecology equipment aid in complete diagnosis of a range of medical problems. Our extensive range of sizes ensures all needs during procedures can easily be met, allowing for better visualisation during an examination leading to improved comfort for the patient.

The smoke evacuation tube enables a clear view.

The high quality construction with secure self retention and durable insulation allow for ease of use for both left & right handed users with reliable longevity.

RB Medical have a large range of speculums available. If you don't see what you're looking for please contact your local area manager who will be more than happy to see what we can do to help.

Our speculums provide a high quality cost effective solution for your surgical needs - if you wish to discuss options available for this product please see our website or contact your local representative.

As a company we have the most exacting standards and comply to the safety and quality standards.

- Complies to regulations ISO 13485 & ISO 9001
- UK engineered to exacting standards
- Complies to the Medical Device Directive (MDD)